

## Cross references: OEM &amp; Alpha Resources Article Numbers

Elementar Article No.	Alpha Resources Article No.	Description
03 001 118	AR-AEB2147	SULPHANILIC ACID, OAS 5 G
03 002 252	AR-02252	O RING SILICONE RUBBER
03 002 253/4	AR-2253/4	REDUCTION TUBE QUARTZ RAPID N
03 002 262	AR-AEB1196	PLATINUM CATALYST 5% Pt 30G
03 002 262	AR-049	N CATALYST, 50 GRAM BOTTLE
03 002 267	AR-02267	O RING VITON
03 002 268	AR-02268	O RING
03 002 395	AR-02395	O-RING GRIPPER ARM
03 002 399	AR-AEE1174	FILTER DIAMETER 18mm PACK OF 1
03 002 607	AR-02607	SHAFT SEAL
03 002 849	AR-01040	POST COMBUSTION TUBE STAINLESS
03 002 850	AR-AEE1176	O RING VITON, 8mm x 2mm
03 651 011	AR-51011	TUBING, PTFE 5 X 3 MM, .5M L
03 653 103	AR-53103	O RING, SILICONE
03 653 118	AR-53118	O RING SILICONE
03 654 505	AR-54505	O RING SILICONE
03 654 603	AR-54603	O RING, SILICONE, PACK OF 10
03 654 615	AR-54615	O RING NITRILE
03 654 627	AR-54627	O RING VITON
03 654 629	AR-54629	O RING, SILICONE
03 654 701	AR-54701	O RING NITRILE
03 654 772	AR-54772	O RING
03 654 814	AR-54814	O RING, SILICONE RUBBER
03 667 990	AR-67990	SEALING RING
03 667 991	AR-67991	SEALING RING
03 674 012	AR-AED5132	ALUMINUM WEIGHING BOATS 12 x 4
03 674 016	AR-AED5133	ALUMINUM WEIGHING BOATS 12 x 6
03 679 900	AR-2174	ALPHASOLVE II, 20-30 MESH 500
03 679 901	AR-171	ALPHADRI MAGNESIUM PERCHLORATE
03 679 905	AR-01029	COPPER OXIDE WIRE, 250 GRAMS,
03 679 907	AR-AEB1008	QUARTZ (SILICA) CHIPS, 50 G
03 679 908	AR-177	QUARTZ WOOL FINE, 4-12 micron
03 679 910	AR-79910	HIGH PURITY CARBON GRANULAR
03 679 914	AR-79914	LEAD CHROMATE - DISCONTINUED
03 679 915	AR-79915	CERIUM DIOXIDE GRANULAR, 100G
03 679 924	AR-189	COPPER STICKS 100 G AMPOULE
03 679 941	AR-79941	O RING SILICONE
03 679 942	AR-1053	ACETANILIDE, 10 G
03 950 005	AR-AEE2291	SPATULA MICRO SPOON FORMED BEN
03 951 619	AR-AED4132	SMOOTH WALL TIN CAPSULES FLAT
03 951 620	AR-AED4131	SMOOTH WALL TIN CAPSULES FLAT
03 960 157	AR-60157	SILICA GEL, 1 KG
03 960 350	AR-AEB2132	NICOTINAMIDE OAS, 1 GRAM
03 960 352	AR-AEB1197	SICAPENT (2.8L) 500G
03 968 006	AR-AEE2240	SEE AR241
03 968 006	AR-241	SILICONE GREASE 150g
03.002.278	AR-2278	METAL CRUCIBLES 100 PACK

### Cross references: OEM & Alpha Resources Article Numbers

Elementar Article No.	Alpha Resources Article No.	Description
05 000 043	AR-AEE1062	O RING SILICONE RUBBER
05 000 043	AR-00043	O RING SILICONE
05 000 044	AR-AEE2518	CAP GL 18
05 000 046	AR-AEE2519	CAP GL25
05 000 049	AR-AEE2517	GL14 SCREW CAP
05 000 086	AR-00086	ALUMINA WOOL FINE, 20G
05 000 087	AR-00087	ASH FINGER CERAMIC RAPID N
05 000 088	AR-50088	OXYGEN LANCE, 140MM
05 000 092	AR-00092	O RING, KALREZ
05 000 094	AR-AEE2505	GLOVES PAIR
05 000 095	AR-00095	O RING VITON
05 000 096	AR-00096	O RING VITON
05 000 140	AR-00140	O RING NITRILE
05 000 247	AR-00247	SEAL, VITON RUBBER
05 000 248	AR-00248	O RING SILICONE
05 000 249	AR-00249	O RING VITON
05 000 308	AR-2167	Smooth Wall Tin Capsule (H=22,
05 000 311	AR-2166	Smooth Wall Tin Capsule (H= 18
05 000 312	AR-AED4120	SMOOTH WALL TIN CAPSULES FLAT
05 000 316	AR-05316	NITROGEN FREE PAPER, 4x4, 500
05 000 341	AR-00341	O RING, VITON RUBBER
05 000 353	AR-50353	SQUARE GASKET FOR ELEMENTAR
05 000 354	AR-00354	OXYGEN LANCE
05 000 359	AR-00359	O RING VITON
05 000 368	AR-00368	O RING
05 000 369	AR-00369	O RING VITON RUBBER
05 000 371	AR-00371	O RING VITON
05 000 372	AR-00372	O RING VITON
05 000 384	AR-00384	QUADRING NITRILE
05 000 385	AR-00385	O RING VITON
05 000 400	AR-AEB2047	BENZOIC ACID OAS 5G
05 000 416	AR-00416	O RING VITON
05 000 425	AR-00425	O RING
05 000 429	AR-2186	TIN FOIL CUPS, 36mm x 36mm
05 000 429	AR-2186-200	TIN FOIL CUPS, 36mm x 36mm
05 000 468	AR-502	BENZOIC ACID POWDER, 99%
05 000 469	AR-00469	CERMIUM DIOXIDE FUSED GRANULAR
05 000 481	AR-0481	Pressed Tin Capsule, (H=11.5mm
05 000 481	AR-059	Smooth Wall Tin Capsule (H=18m
05 000 568	AR-AEE1175	O RING 16mm X 3mm
05 000 658	AR-AEC3013	OXYGEN LANCE 170MM
05 000 658	AR-00658	OXYGEN LANCE, 170 MM
05 000 699	AR-189	COPPER STICKS 100 G AMPOULE
05 000 814	AR-00814	ALUMINUM CAPSULES SMOOTH WALL,
05 000 815	AR-00815	ALUMINUM CAPSULES SMOOTH WALL
05 000 816	AR-00816	SILVER CAPSULES, SMOOTH WALL

### Cross references: OEM & Alpha Resources Article Numbers

Elementar Article No.	Alpha Resources Article No.	Description
05 000 880	AR-060	MANGANESE DIOXIDE
05 000 900	AR-00900	Quartz Chips (Turnings) , 2 TO
05 001 013	AR-338	ALPHAGEL CAPSULES, 7 MM OD, 40
05 001 039	AR-01039	COPPER OXIDE WIRE, FINE WIRES,
05 001 053	AR-01053	O RING FOR SPRING FLANGE
05 001 059	AR-AEB1141	MAGNESIUM PERCHLORATE SELF IND
05 001 067	AR-01067	GLUTAMIC ACID, 100 GRAMS
05 001 070	AR-AEB1002	COPPER OXIDE, FINE WIRES, .5 X
05 001 070	AR-01070	COBALTOUS/IC OXIDE GRANULAR, .
05 001 081	AR-2019	WHEAT FLOUR C= 38.86% H= 5.45%
05 001 169	AR-AEB2148	PHENYLALANINE OAS 25G
05 001 432	AR-AEE2520	BRASH BRUSH
05 001 715	AR-01715	COTTON WOOL, 50 GRAMS
05 002 321	AR-AEC3028	OXYGEN LANCE 90mm
05 003 394	AR-ATD1008	PRESSED TIN CAPSULES 8 X 5, 25
05 003 400	AR-059	Smooth Wall Tin Capsule (H=18m
11 001 170/4	AR-AEC3018	ASH FINGER
11 001 171/4	AR-AEC3016	SUPPORT TUBE
11 001 171/4	AR-11714	SUPPORT TUBE, L=133MM
11.00-1150/4	AR-11504	COMBUSTION/REDUCTION TUBE
11.00-1170/4	AR-1170/4	ASH FINGER FOR ELEMENTAR
11.00-1308/4	AR-AEC3010	REDUCTION TUBE VARIO-EL
11.00-1308/4	AR-13084	REDUCTION TUBE
11.00-1317/4	AR-AEC3021	SEE AR1317/4
11.00-1317/4	AR-1317/4	PROTECTIVE TUBE FOR ASH FINGER
11.00-1317/4	AR-13174	REPLACED BY AR1317/4
11.00-1319/4	AR-AEC3068	ASH CRUCIBLE, QUARTZ, 30MM
11.00-1331/4	AR-AEC3022	CONNECTING TUBE TRANSPARENT SI
11.00-1331/4	AR-13314	CONNECTING TUBE TRASPARENT SIL
11.00-1371	AR-01371	SILVER WEIGHING BOATS, 12X4X4
11.00/1172/4	AR-1172/4	QUARTZ SUPPORT TUBE
11.02-0008	AR-2181	TUNGSTEN TRIOXIDE 60 GRAM
11.02-0017	AR-20017	TUNGSTEN OXIDE FINE POWDER
11.02-0024	AR-20024	TUNGSTEN OXIDE POWDER, 90G
11.02-1036/4	AR-11364	ASH FINGER
11.02-1047	AR-AED5029	TIN WEIGHING BOATS 6 x 4 x 4mm
11.22-0003/3	AR-AEC3056	COMBUSTION TUBE VARIO TRACE S
11.45-1005/4	AR-AEC3062	ASH FINGER, STAINLESS STEEL, 9
11.45-1012	AR-AEC3014	OXYGEN LANCE 190mm
11.45-1012/4	AR-10124	OXYGEN LANCE, 190 MM
11001317/4	AR-AEC3059	PROTECTION TUBE VARIO EL
12 001 008	AR-929	CRUCIBLE TONGS FOR AR3818
12 010 032	AR-10032	TIN FOIL SQUARES, 36X36 MM, PA
12.00-0040	AR-AEB1205	TUNGSTEN METAL, 500 G, FOR USE
12.00-0040	AR-02530	TUNGSTEN METAL, 500G, FOR USE
12.00-0040	AR-02530B	SEE AR02530 or AR02530-2500

### Cross references: OEM & Alpha Resources Article Numbers

Elementar Article No.	Alpha Resources Article No.	Description
12.00-1004/4	AR-01004/4	WIRE SCREEN, 22 MM, FOR ELEMEN
12.00-1005/4	AR-AEC3035	WIRE MESH RONDE 18mm
12.00-1005/4	AR-10054	WIRE MESH ROUND, 18 MM
12.00-1006/4	AR-01020	COMBUSTION TUBE STAINLESS STEE
12.00-1009	AR-01009	TIN BOATS, 9X10X20MM, 120/PK
12.00-1087-4	AR-10874	ASH FINGER STAINLESS STEEL
12.00-1087/4	AR-AEC3031	ASH FINGER STAINLESS STEEL
12.00.0040	AR-266	ALPHACEL TUNGSTEN -12+20
12.01-0033	AR-12014	TIN FOIL SQUARES, 50X50 MM, 10
12.01-0033	AR-2397	TIN FOIL CUPS, 51mm x 51mm
12.01-0033	AR-2397-400	TIN FOIL CUPS, 51mm x 51mm
12.01-0045	AR-02530X	TUNGSTEN KIT COMPLETE WITH 500
12.01-1008/4	AR-02530S	SPACERS FOR TUNGSTEN, PACK OF
12.01-1008/4	AR-02530S-200	SPACERS FOR ELEMENTAR, FOR USE
12.01-1008/4	AR-02530S-500	SPACERS FOR ELEMENTAR, FOR USE
12.01-1012/4	AR-10124	OXYGEN LANCE, 190 MM
12.01-1013/4	AR-02510	QUARTZ CATALYST TUBE, FOR ELEM
12.01-1032/4	AR-AEC3048	DISTANCE RING
12.01-1032/4	AR-11032/4	DISTANCE RING, FOR ELEMENTAR
12.01-1036/4	AR-AEC3041	COMBUSTION TUBE S/S
12.01-1036/4	AR-10364	COMBUSTION AND POST COMBUSTION
15.00-0062	AR-AEB2048	SULPHANILAMIDE OAS 10G
15.00-0080	AR-AEC1143	DRYING TUBE
15.00-1189	AR-01189	ASH FINGER HIGH TEMPERATURE PY
15.00-1319	AR-1319	ASH CRUCIBLE CERAMIC, NO SLOTS
15.40-1002	AR-AEC3057	SUPPORT ROD REDUCTION TUBE 110
16.00-1012	AR-AEC3026	DRYING TUBE 348mm
16.00-1037	AR-AEC3047	DRYING TUBE SHORT
16.00-1104	AR-AEC3065	ASH CRUCIBLE, STAINLESS STEEL
16.02-0036	AR-AEC3069	COMBUSTION TUBE
22 127 245/3	AR-72453	REDUCTION TUBE, RAPID CHN
22 127 313/4	AR-73134	RECEPTACLE INSERT, CHN RAPID
22 131 365	AR-AEB1231	SEE AR22365
22 133 210/3	AR-32103	PYROLISI TUBE, VARIO-EL
22 133 211/4	AR-32114	PYROLYSIS TUBE INSERT, VARIO-E
22 133 212	AR-33212	SILVER WEIGHING BOAT 12X6X6MM,
22 133 213	AR-33213	SILVER BOATS, 15X8X8MM, 70/PAC
22 133 214	AR-33214	SILVER BOATS, 20X10X9MM, 42/PK
22 137 418	AR-37418	TIN WEIGHING BOAT, 12X4X4MM, 1
22 137 419	AR-37419	TIN WEIGHING BOATS 12X6X6MM, 5
22 137 420	AR-37420	TIN BOATS 8X8X15MM, 250/PK
22 156 202/3	AR-AEC3004	COMBUSTION TUBE RAPID CHN
22 160 200/4	AR-02004	DRYING GUARD TUBE
22 160 204/4	AR-AEC3008	DRYING TUBE
22 160 204/4	AR-AEC3066	DRYING TUBE, STRAIGHT
22 160 301/4	AR-AEC3037	U TUBE RAPID N COMPLETE WITH C

Cross references: OEM & Alpha Resources Article Numbers

Elementar Article No.	Alpha Resources Article No.	Description
22 160 301/4	AR-01032	U TUBE, 195 MM LENGTH
23.00-1020	AR-AEC7017	GLASSY CARBON TUBE PYROCUBE
23.00-1021	AR-01021	CERAMIC TUBE, PYROCUBE
23.00-1148	AR-AEC3052	SUPPORT ROD REDUCTION TUBE 175
23.00-1149	AR-01149	REACTION TUBE, PYROCUBE
23.00-1150	AR-01150	COMBUSTION TUBE, PYROCUBE
23.00-1189	AR-AEC3061	ASH CRUCIBLE, GRAPHITE
25.00-0017	AR-01007	QUARTZ COMBUSTION TUBE, FOR EL
25.00-0029	AR-00029	GRIPPER ARM ASSEMBLY, WITH ARO
25.00-1028	AR-01028	COMBUSTION TUBE, QUARTZ, VARIO
25.00-1030	AR-01030	COMBUSTION TUBE VARIOMAX, STAI
25.00-1048	AR-01048	POST COMBUSTION TUBE
25.20-0004	AR-AEE2522	GRIPPER ARM WITH O RINGS FOR V
25.20-1003	AR-01003	CONNECTING TUBE TRANSPARENT SI
29.00-0092	AR-AEB1407	Zinc Granules
29.00-1078	AR-AEC3067	COMBUSTION TUBE
33.00-2125	AR-02125	REACTOR TUBE, LIQUI-TOC II
34 001 339/4	AR-2216/4	U TUBE, 180MM LENGTH ~;-
34.00-1069/4	AR-AEC3100	COMBUSTION TUBE HIGH TOC II
35.00-2125	AR-AEC3101	REACTOR TUBE liquiTOC
38.00-1290	AR-AEC3054	ASH CRUCIBLE, QUARTZ
50 008 449/4	AR-08449/4	WIRE SCREEN, 26.8 MM
50 008 467	AR-01008	CORUNDUM BALLS , 40g PER BOTTL
50 008 467	AR-01008-200	CORUNDUM BALLS, 355 cc BY VOLU
50 009 375/4	AR-AEC3045	FILTER TUBE
88 895 805	AR-AEB2138	ASPARTIC ACID OAS 100G
88 895 805	AR-01011	L-ASPARTIC ACID, 98%+, OAS, 10
905.200.970.001	AR-2174	ALPHASOLVE II, 20-30 MESH 500
905.200.970.001 AN	AR-2174	ALPHASOLVE II, 20-30 MESH 500