

*Cross references: OEM & Alpha Resources Article Numbers*

VELP	Alpha Resources	Description
10003529	AR-140020	TOP REACTOR GASKET
10005138	AR-AEC1008	SCRUBBER TUBE WATER THREADED B
10005506	AR-140019	LOWER REACTOR GASKET
A00000148	AR-AEB3015	CHROMOSORB W 30-60 MESH ACID W
A00000148	AR-AEB3017	CHROMOSORB W 30-60 MESH ACID W
A00000149	AR-2092	EDTA, 50 G BOTTLE
A00000153	AR-2186	TIN FOIL CUPS, 36mm x 36mm
A00000153	AR-2186-200	TIN FOIL CUPS, 36mm x 36mm
A00000154	AR-177	QUARTZ WOOL FINE, 4-12 micron
A00000155	AR-AEB1091M	COPPER WIRES 3X100G PACK FINE
A00000155	AR-AEB1176	COPPER WIRES FINE WIRES REDUCE
A00000157	AR-AEB1002	COPPER OXIDE, FINE WIRES, .5 X
A00000157	AR-AEB1125	COPPER OXIDE FINE WIRES 4 X 0.
A00000158	AR-140026	PRE PACKED COMBUSTION TUBE FOR
A00000159	AR-140011	ALPHACAT HT CATALYST FOR DUMAT
A00000160	AR-140014	ALPHACAT LT CATALYST FOR DUMAT
A00000161	AR-AEC1131	TURBO N/DUMATHERM QUARTZ INSERT
A00000162	AR-0000162	QUARTZ RECTION TUBE VELP NDA70
A00000162	AR-140018	QUARTZ REACTION TUBES DUMATHER
A00000171	AR-AEB1197	SICAPENT (2.8L) 500G
A00000198	AR-5000	POROUS RETICULATED CRUCIBLE,
A00000225	AR-171	ALPHADRI MAGNESIUM PERCHLORATE
A00000226	AR-140226	PRE-PACKED COPPER REDUCTION TU
A00000235	AR-2028	RICE FLOUR C= 45.38% H= 6.69%
A00000236	AR-241	SILICONE GREASE 150g
A00000240	AR-140007	COPPER STICK FOR DUMATHERM/COS
A00000240	AR-140246	HIGH REDUCTION COPPER STICKS -