

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
002-301	AR-3818	CERAMIC CRUCIBLE, BAGGED, BOX
002-302	AR-3818	CERAMIC CRUCIBLE, BAGGED, BOX
002-308	AR-3818F	CERAMIC CRUCIBLE, BOX OF 1000,
002-309	AR-3818F	CERAMIC CRUCIBLE, BOX OF 1000,
145-601	AR-5601	WATER TEMP REGULATOR ~;~
156-002	AR-002	O RING UPPER ELECTRODE
156-003	AR-156	O RING
156-004	AR-005	BACKUP RING, TEFLON ~;~
163-012	AR-63012	TUBING INSERT
170-001	AR-0001	SPRING, LOWER ELECTRODE ~;~
172-505	AR-2505	TEFLON TUBING 1/4" OD X 3/16" I
175-573	AR-5573	DISCONTINUED - NO REPLACEMENT
175-574	AR-5574	OPTICAL SWITCH ASSEMBLY
175-590	AR-5590	PUMP REPAIR KIT
175-631	AR-5631	MOTOR BRUSH FOR 774-571 VACUUM
175-632	AR-5632	Flow Controller "O" Ring - ALL
175-659	AR-5659	PAR-APER DEFLECTOR, V55147 ~;~
175-673	AR-5673	NO LONGER AVAILABLE!!!
175-674	AR-5674	TEAR OFF BLADE SPRING
183-001	AR-2174	ALPHASOLVE II, 20-30 MESH 500
188-012	AR-293	BRUSH, CLEANING, SMALL
190-142	AR-142	SS Electrode Plate Screws, pkg
190-322	AR-AR-AEE2248	SCREW FOR FRONT U BLOCK
190-739	AR-739	NYLON SPACER SCREW
190-742	AR-AR-AEE2241	3/8 SOCKET CAR-AP SCREW
194-060	AR-4060	HEX JAM NUT
194-150	AR-4060	HEX JAM NUT
204-001	AR-001	THERMOREGULATOR/TC30/TN14
501-001	AR-1701	0.57% SULFUR COAL STANDARD. 50
501-005	AR-1704	1.09% SULFUR IN COAL STANDARD,
501-008	AR-008	ALPHACEL II TUNGSTEN/TIN, -20/
501-009	AR-1009	SPRING CLIP FOR DUST TRAR-AP ~;~
501-010	AR-010	DUST TRAR-AP/COMPLETE
501-011	AR-011	DUST TRAR-AP INSIDES ~;~
501-012	AR-011B	DUST RESERVOIR ~;~
501-013	AR-013	DUST TRAR-AP FILTER CLOTH, PKG OF
501-014	AR-1705	1.48% SULFUR IN COAL
501-015	AR-1703	0.89% SULFUR COAL STANDARD, 50
501-020	AR-1708	3.03% SULFUR IN COAL STANDARD
501-022	AR-1709	3.79% SULFUR IN COAL STANDARD
501-024	AR-322	CAST IRON STANDARD, C= 3.38% S
501-027	AR-204	RUBBER SQUEEZE BULB
501-031	AR-031	INDUCTION COIL L TYPE
501-034	AR-1034	CALCIUM CARBONATE STANDARD, 12
501-034	AR-1034-100	CALCIUM CARBONATE STANDARD, 12
501-045	AR-1045	ZIRCON CRUCIBLE, CHN
501-045-125	AR-1045	ZIRCON CRUCIBLE, CHN

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
501-046	AR-2021	SUCROSE STANDARD, 50 GRAM
501-050	AR-AEB2099	NICOTINIC ACID OAS 15G
501-053	AR-AEB2000	ACETANILIDE OAS 1g
501-053	AR-1053	ACETANILIDE, 10 G
501-057	AR-079G	ALPHANAL WETTING AGENT, 1gal.
501-059	AR-059	Smooth Wall Tin CAR-APsule (H=18m
501-060	AR-060	MANGANESE DIOXIDE
501-061	AR-1061	ARROWROOT STARCH, 227g
501-062	AR-061	BOAT PULLER
501-067	AR-067	ROTOMETER BALL ~;-~
501-073	AR-073	GRAR-APHITE POWDER, 100 G
501-073	AR-539	DISSONTINUED
501-076	AR-076	Tin Metal Accelerator 20 Mesh
501-077	AR-077	IRON CHIP ACCELERATOR, 908 G
501-078	AR-078	IRON POWDER, 1 POUND CONTAINER
501-079	AR-079	ALPHANAL WETTING AGENT 4oz
501-081	AR-081	GLASS WOOL, 1 POUND CONTAINER
501-082	AR-022	CLEANING BRUSH
501-083	AR-304	DUCTILE IRON C= 4.21% S= 0.021
501-084	AR-084	TYGON TUBING 1/4" ID, 3/8" OD
501-085	AR-086	TYGON TUBING 5/16"" ID
501-087	AR-087	TYGON TUBING 3/16"" ID
501-089	AR-507	ITEM DISCONTINUED - NO REPLACE
501-095	AR-2814	RESIDUAL FUEL OIL STANDARD 1.0
501-099	AR-2837	LUBE OIL, 4.00%S , 100ml
501-101	AR-800	ARROWROOT STARCH SOLUTION, 500
501-105	AR-301	CAST IRON C=1.99% S=0.0162%,
501-110	AR-110	SILICONE TUBING, 9.5 MM ID X 3
501-110-009	AR-110	SILICONE TUBING, 9.5 MM ID X 3
501-110-019	AR-110-19	SILICON TUBING, 3/8" ID, 19" L
501-112	AR-112	1/4"" I.D. NEOPRENE TUBING
501-123	AR-307	TUNGSTEN CARBIDE C= 6.13% S= 0
501-130	AR-1710	4.25% SULFUR IN COAL STANDARD
501-130	AR-1714	1.71% SULFUR IN COAL STANDARD,
501-144	AR-144	MAGNESIUM OXIDE, 1 POUND
501-147	AR-148	COPPER PIN STANDARD, HIGH OXYG
501-148	AR-149	COPPER PIN STANDARD, HIGH OXYG
501-165	AR-1165	DRIERITE INDICATING 10-20 mesh
501-166	AR-166	REPLACED BY AR 2176
501-170	AR-170	RARE EARTH COPPER OXIDE, 50 GR
501-171	AR-171	ALPHADRI MAGNESIUM PERCHLORATE
501-171-HAZ	AR-171	ALPHADRI MAGNESIUM PERCHLORATE
501-179	AR-079P	ALPHANAL WETTING AGENT, 1pint
501-184	AR-1711	5.13% SULFUR IN COAL STANDARD
501-197	AR-197	REPLACED BY AR 775
501-233	AR-233	CERAMIC ADHESIVE ALUMINA, REP
501-234	AR-234	OXYGEN SPECIMEN FILE

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
501-241	AR-241	SILICONE GREASE 150g
501-263	AR-263	COPPER ACCELERATOR, 3 POUNDS
501-264	AR-264	COPPER ACCELERATOR, 227 GRAMS
501-265	AR-1265	CRUCIBLE HOLDER TRAY
501-280	AR-634	""O"" RING (DUST TRAR-AP)
501-291	AR-291	OXYGEN REGULATOR
501-369	AR-1727	PROX COAL STANDARD, 50 GRAMS
501-426	AR-426	COMBUSTION CONTROL, SOLID, I P
501-427	AR-427	COMBUSTION CONTROL FOR LIQUID
501-428	AR-1721	PROX COAL STANDARD, 50 GRAM BO
501-429	AR-1723	PROX COAL STANDARD, 50 GRAM BO
501-430	AR-1724	PROX COAL STANDARD, 50 GRAM BO
501-431	AR-1726	PROX COAL STANDARD, 50 GRAMS
501-432	AR-1728	PROX COAL STANDARD, 50 GRAMS
501-433	AR-1730	PROX COAL STANDARD, 50 GRAM BO
501-438	AR-734	PROX MET COKE, 50 GRAM
501-439	AR-000	PARAFFIN OIL C,H,N, BTU
501-441	AR-2021	SUCROSE STANDARD, 50 GRAM
501-457	AR-400	TIN COATED COPPER, 3 POUND CON
501-458	AR-2822	#2 DIESEL OIL 0.305% S, 100ml
501-501	AR-870	C/S RING STANDARD C= 0.027% S
501-502	AR-871	C/S RING STANDARD C= 0.0363% S
501-503	AR-872	C/S RING STANDARD C= 0.159% S=
501-504	AR-873	C/S RING STANDARD C=0.428% S=0
501-505	AR-874	C/S RING STANDARD C= 0.658% S
501-506	AR-875	C/S RING STANDARD C= 0.799% S=
501-529	AR-556	HYDROGEN IN STEEL PIN STANDARD
501-539	AR-930	ACTIVATED CHARCOAL - ITEM DISC
501-551	AR-660	O/N IN STEEL PIN STANDARD O=0.
501-551	AR-667	O/N IN STEEL PIN STANDARD O= 0
501-551	AR-670	O/N IN STEEL PIN STANDARD O= 0
501-552	AR-656	O/N IN STEEL PIN STANDARD O=0.
501-552	AR-662	O/N IN STEEL PIN STANDARD O= 0
501-553	AR-669	SOLD OUT & DISCONTINUED. NO R
501-561	AR-2016	SOY BEAN MEAL N= 8.40% S= 0.46
501-563	AR-2025	CORN MEAL STANDARD. 30g PER B
501-571	AR-571	COPPER CAR-APSULE - bottle of 100
501-587	AR-912	PLATINIZED SILICA
501-587	AR-912B	PLATINIZED SILICA, BULK 500G U
501-592	AR-062	MANGANESE DIOXIDE
501-598	AR-598	SHORT NIBBLED NICKEL LOW OXYGE
501-598	AR-599	LONG NIBBLED NICKEL LOW OXYGEN
501-608	AR-1608	QUARTZ WOOL WADDING
501-609	AR-615	FURNACE REAGENT, 100g PER BOTT
501-614	AR-1614	SPATULA, RUBBER GRIP HANDLE, F
501-621	AR-621	COPPER TURNINGS, DEGASSED, 200
501-632	AR-632	MOLE SIEVE TRAR-AP 14 G

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
501-636	AR-636	VANADIUM PENTOXIDE, 250 G BOTT
501-640	AR-366	COPPER ACCELERATOR, 25 POUND
501-641	AR-643	COPPER METAL SCRUBBER
501-644	AR-644	O/N IN STEEL PIN STANDARD O= 3
501-644	AR-645	O/N IN STEEL PIN STANDARD O= 0
501-645	AR-1651	O/N IN STEEL PIN STANDARD O= 0
501-647	AR-1655	SOLD OUT NO REPLACEMENT
501-648	AR-2776	ULTIMATE COAL STANDARD, 50 GRA
501-656	AR-621	COPPER TURNINGS, DEGASSED, 200
501-673	AR-673	WROUGHT IRON CHIP STANDARD C=0
501-674	AR-881	C/S PIN STANDARD C= 0.0161% S
501-675	AR-882	C/S PIN STANDARD C= 0.0532% S=
501-675	AR-891	C/S PIN STANDARD C= 0.0236% S=
501-676	AR-890	STEEL PIN STANDARD CARBON= 0.0
501-677	AR-895	C/S PIN STANDARD C= 0.275% S=
501-678	AR-884	C/S PIN STANDARD C= 0.443% S=
501-678	AR-886	C/S PIN STANDARD C= 0.724% S=
501-679	AR-885	C/S PIN STANDARD C= 0.821% S=
501-679	AR-889	C/S PIN STANDARD C= 0.931% S=
501-697	AR-612	PLATINIZED CHARCOAL ~;~
501-951	AR-1951	NICKEL CALIBRATION WIRE
501-955	AR-1722	PROX COAL STANDARD, 50g
501-957	AR-1729	PROX COAL STANDARD, 50 GRAM BO
501-961	AR-304	DUCTILE IRON C= 4.21% S= 0.021
501-964	AR-1712	5.79% SULFUR IN COAL STANDARD
501-964	AR-1715	7.06% SULFUR IN COAL STANDARD,
501-967	AR-2851	EXTREME PRESSURE OIL , 100ml
501-968	AR-2852	EXTREME PRESSURE OIL 2.95% S,
501-969	AR-2853	EXTREME PRESSURE OIL , 100ml
502-005	AR-1205	PIPETTE BRUSH
502-007	AR-2189	TUBE CLEANER, 8MM X 4FT, FP 42
502-008	AR-570	TIN PLUGS, 150 PER BOTTLE
502-016	AR-655	SOLD OUT - REPLACEMENT EXPECTE
502-016	AR-663	O/N IN STEEL PIN STANDARD O= 0
502-020	AR-599	LONG NIBBLED NICKEL LOW OXYGEN
502-021	AR-548	DISCONTINUED ~;~
502-023	AR-473	UTILITY FUNNEL FOR GLASS TUBES
502-024	AR-651	HON IN TITANIUM STANDARD H=0.0
502-026	AR-2815	RESIDUAL FUEL OIL 2.03% SULFUR
502-026	AR-2818	RESIDUAL FUEL OIL STANDARD 1.5
502-028	AR-1617	TWEEZERS, STAILNESS STEEL, STR
502-032	AR-078	IRON POWDER, 1 POUND CONTAINER
502-040	AR-2166	Smooth Wall Tin CAR-APsule (H= 18
502-042	AR-025	QUARTZ CRUCIBLE, CHN
502-049	AR-049	N CATALYST, 50 GRAM BOTTLE
502-052	AR-323	PIG IRON STANDARD, C= 3.92% S
502-055	AR-AEB2141	ORCHARD LEAVES STANDARD 20G

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
502-060	AR-556	HYDROGEN IN STEEL PIN STANDARD
502-061	AR-555	HYDROGEN IN STEEL PIN STANDARD
502-062	AR-AEB2180	SOIL STANDARD (SANDY) 30G
502-062	AR-4018	C/S IN SOIL STANDARD. C= 1.26%
502-062	AR-4020	CARBON IN SYNTHETIC SOIL STAND
502-062	AR-4029	CARBON IN SYNTHETIC SOIL STAND
502-062	AR-4105	LOSS ON IGNITION STANDARD 5% L
502-062	AR-4106	LOSS ON IGNITION STANDARD 10%
502-062	AR-4107	LOSS ON IGNITION STANDARD 20%
502-062	AR-4108	LOSS ON IGNITION STANDARD 30%
502-071	AR-2864	LUBE OIL 29.67% S , 100ml
502-072	AR-655	SOLD OUT - REPLACEMENT EXPECTE
502-076	AR-571	COPPER CAR-APSULE - bottle of 100
502-077	AR-570	TIN PLUGS, 150 PER BOTTLE
502-083	AR-100	RESIDUAL OIL, C,H,N, BTU 4oz
502-085	AR-3085	ZINC SULFIDE STANDARD 32%, 50
502-091	AR-091	CALCIUM OXALATE 50g
502-092	AR-2092	EDTA, 50 G BOTTLE
502-092-250	AR-2092-250	EDTA, 250G BOTTLE
502-122	AR-2122	FIRING FUSE WIRE, 2" LONG
502-134	AR-2134	BENZOIC ACID TABLET,1.45G, PAC
502-134	AR-2134C	BENZOIC ACID TABLET, 1.45G/TAB
502-134	AR-2134V	BENZOIC ACID TABLET,1.45G, PAC
502-154	AR- 648	SEE AR648
502-154	AR-648	H/O/N IN TITANIUM; H= 0.0155%
502-156	AR-2156	FLUOR HIBIT, 1 POUND CONTAINER
502-167	AR-2167	Smooth Wall Tin CAR-APsule (H=22,
502-173	AR-173	ALPHACEL II HP, 750G BOTTLE
502-174	AR-167	REPLACED BY AR2174
502-174	AR-2174	ALPHASOLVE II, 20-30 MESH 500
502-176	AR-169	ALPHASOLVE 8-20 MESH, 250 G UN
502-176	AR-2176	ALPHASOLVE II, 8-20 MESH 500 G
502-177	AR-177	QUARTZ WOOL FINE, 4-12 micron
502-181	AR-2181	TUNGSTEN TRIOXIDE 60 GRAM
502-183	AR-2183	NICKEL CAR-APSULE, 6.5 X 10 MM, F
502-184	AR-502	BENZOIC ACID POWDER, 99%
502-186	AR-2186	TIN FOIL CUPS, 36mm x 36mm
502-186	AR-2186-200	TIN FOIL CUPS, 36mm x 36mm
502-186-100	AR-2186	TIN FOIL CUPS, 36mm x 36mm
502-186-100	AR-2186-200	TIN FOIL CUPS, 36mm x 36mm
502-186-1000	AR-2186-1000	Tin Foil Cups, 36 x 36mm, 1000
502-186-200	AR-2186-200	TIN FOIL CUPS, 36mm x 36mm
502-187	AR-2187	Smooth Wall Tin CAR-APsule (H=9mm
502-188	AR-2188	ALUMINUM OXIDE PELLETS
502-189	AR-189	COPPER STICKS 100 G AMPOULE
502-190	AR-AEB1002	COPPER OXIDE, FINE WIRES, .5 X
502-202	AR-1701	0.57% SULFUR COAL STANDARD. 50

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
502-203	AR-AEB2031	STEARIC ACID STANDARD, OAS
502-204	AR-AEB2034	3,5-DINITROBENZOIC ACID, OAS
502-205	AR-AEB2030	CAFFEINE STANDARD, OAS
502-206	AR-ASD2000	SILVER CAR-APSULES, PRESSED, 3.2
502-206	AR-ASD2003	AG CAR-APSULES, 5MM X 3.5MM PACK
502-207	AR-AEB2035	CYSTINE, OAS
502-208	AR-AEB2036	SULPHANILAMIDE, OAS,1 GRAM BTL
502-209	AR-2209	SULFAMETHAZINE, 5 G
502-211	AR-211	GLYCINE
502-213	AR-2213	LEAK DETECTION SOLUTION, 32 ou
502-225	AR-AEB1036	AG TUNGSTATE/MG OXIDE
502-226	AR-AEB1009	AG VANADATE REAGENT, 10g PER B
502-227	AR-2227	PLEASE SEE ATD-1000 ~;~
502-227	AR-ATD1000	TIN CAR-APSULES, PRESSED, 4 X 3.2
502-231	AR-673	WROUGHT IRON CHIP STANDARD C=0
502-258	AR-621	COPPER TURNINGS, DEGASSED, 200
502-266	AR-205	FIRING FUSE WIRE, 300 DETERMIN
502-272	AR-2017	CORN GLUTEN C= 45.27% H= 5.98%
502-273	AR-AEB2018	4TRIFLUOROMETHYLBENZAMIDE, OAS
502-273	AR-2018	ALFALFA N= 3.52% S= 0.27%, 30
502-274	AR-2019	WHEAT FLOUR C= 38.86% H= 5.45%
502-275	AR-2020	RYE FLOUR C= 39.50% H= 5.50% N
502-276	AR-2026	OAT MEAL C= 47.00% H= 6.67% N=
502-277	AR-2027	BARLEY FLOUR, 30g PER BOTTLE
502-278	AR-2028	RICE FLOUR C= 45.38% H= 6.69%
502-282	AR-189	COPPER STICKS 100 G AMPOULE
502-283	AR-AEE2105	ABSORBENT DISCS 3.5mm PACK OF
502-294	AR-008	ALPHACEL II TUNGSTEN/TIN, -20/
502-294	AR-173	ALPHACEL II HP, 750G BOTTLE
502-295	AR-621	COPPER TURNINGS, DEGASSED, 200
502-297	AR-008B	ALPHACEL II TUNGSTEN/TIN, -20/
502-298	AR-298	SULFAMETHAZINE STANDARD, 25 G
502-301	AR-2301	TUNGSTEN OXIDE, 125g PER BOTTL
502-304	AR-2304	COPPER STICKS, RESEALABLE BOTT
502-304-500	AR-2304-500	COPPER STICKS, RESEALABLE CONT
502-308	AR-AEB2178	MEDIUM ORGANIC CONTENT SOIL ST
502-308	AR-4016	C/S IN SOIL STANDARD. C= 2.24%
502-308	AR-4021	CARBON IN SYNTHETIC SOIL STAND
502-309	AR-AEB2176	SOIL STANDARD (PEATY) 30G, C=1
502-309	AR-4012	C/S IN LIMESTONE STANDARD C= 1
502-310	AR-2310	STEEL WOOL FILTERING MATERIAL
502-311	AR-2311	FURNACE REAGENT 1000 CALCIUM T
502-314	AR-177	QUARTZ WOOL FINE, 4-12 micron
502-318	AR-4017	C/S IN SOIL STANDARD. C= 0.50%
502-319	AR-4005	ORE STANDARD C= 1.42% S= 1.33%
502-320	AR-4006	ORE STANDARD C= 3.75% S= 3.82
502-321	AR-321	ALPHACAT 250g PER BOTTLE

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
502-335	AR-2335	DISCONTINUED
502-338	AR-338	ALPHAGEL CAR-APSULES, 7 MM OD, 40
502-343	AR-2343	NICKEL BOAT LINER, FOR LIQUID
502-344	AR-2344	NI BASKETS LOW OXYGEN, 1g, 10
502-345	AR-2345	NICKEL BASKETS 1.5G BASKETS, 1
502-351	AR-351	HALOGEN SCRUBBER, 50 G
502-359	AR-2359	MAGNESIUM OXIDE FURNACE REAGEN
502-372	AR-4007	SOLD OUT - Replacement target
502-374	AR-2374	PARTICULATE FILTER, PKG OF 50
502-382	AR-339	ALPHAGEL CAR-APSULES, 8 MM OD, 40
502-383	AR-1700	0.35% SULFUR IN COAL STANDARD,
502-384	AR-1702	0.73% SULFUR COAL STANDARD
502-385	AR-1704	1.09% SULFUR IN COAL STANDARD,
502-386	AR-1705	1.48% SULFUR IN COAL
502-387	AR-1708	3.03% SULFUR IN COAL STANDARD
502-388	AR-1711	5.13% SULFUR IN COAL STANDARD
502-389	AR-2892	RESIDUAL FUEL OIL STANDARD 0.2
502-390	AR-2895	RESIDUAL FUEL OIL STANDARD 0.4
502-391	AR-2814	RESIDUAL FUEL OIL STANDARD 1.0
502-392	AR-2815	RESIDUAL FUEL OIL 2.03% SULFUR
502-393	AR-2816	RESIDUAL FUEL OIL 3.02% SULFUR
502-394	AR-2899	4.50%S IN RESIDUAL OIL, 100ml
502-397	AR-2397	TIN FOIL CUPS, 51mm x 51mm
502-397-400	AR-2397-400	TIN FOIL CUPS, 51mm x 51mm
502-417	AR-2838	LUBE OIL, .01%S , 100ml
502-418	AR-2932	LUBE OIL, 0.05 %S , 100ml
502-419	AR-2934	LUBE OIL, 0.10 %S , 100ml
502-420	AR-2833	LUBE OIL, .50%S , 100ml
502-421	AR-2835	LUBE OIL, 1.00%S , 100ml
502-422	AR-2836	LUBE OIL, 2.00%S, 100 ML BOTTL
502-423	AR-2837	LUBE OIL, 4.00%S , 100ml
502-424	AR-2871	ULTRA LOW #2 DIESEL OIL, NOMIN
502-425	AR-2873	ULTRA LOW #2 DIESEL OIL, NOMIN
502-426	AR-2875	ULTRA LOW #2 DIESEL OIL, NOMIN
502-427	AR-2877	ULTRA LOW #2 DIESEL OIL, NOMIN
502-428	AR-2823	#2 DIESEL OIL, 0 .53% S, 100ml
502-429	AR-2825	#2 DIESEL OIL 1.01% S, 100 ML
502-430	AR-2827	#2 DIESEL OIL 2.00% S, 100ml
502-431	AR-2829	#2 DIESEL OIL 2.95% S, 100ml
502-433	AR-1700	0.35% SULFUR IN COAL STANDARD,
502-433	AR-1701	0.57% SULFUR COAL STANDARD. 50
502-434	AR-1702	0.73% SULFUR COAL STANDARD
502-435	AR-1704	1.09% SULFUR IN COAL STANDARD,
502-436	AR-1705	1.48% SULFUR IN COAL
502-436	AR-1706	2.07% SULFUR COAL STANDARD
502-437	AR-1708	3.03% SULFUR IN COAL STANDARD
502-438	AR-1710	4.25% SULFUR IN COAL STANDARD

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
502-438	AR-1711	5.13% SULFUR IN COAL STANDARD
502-438	AR-1714	1.71% SULFUR IN COAL STANDARD,
502-450	AR-450	HEATED PARTICULATE FILTER, 33M
502-473	AR-473	UTILITY FUNNEL FOR GLASS TUBES
502-489	AR-2489	MICRO POROUS CRUCIBLE FOR CH
502-493	AR-2156	FLUOR HIBIT, 1 POUND CONTAINER
502-634	AR-2029	HIGH PURITY GRAR-APHITE
502-642	AR-AEB2148	PHENYLALANINE OAS 25G
502-654	AR-AEB2049	2,5-(Bis(5-tert-butyl-2-benzo-
502-656	AR-621	COPPER TURNINGS, DEGASSED, 200
502-670	AR-1701	0.57% SULFUR COAL STANDARD. 50
502-673	AR-1708	3.03% SULFUR IN COAL STANDARD
502-674	AR-1710	4.25% SULFUR IN COAL STANDARD
502-674	AR-1714	1.71% SULFUR IN COAL STANDARD,
502-705	AR-189	COPPER STICKS 100 G AMPOULE
502-810	AR-810	ALPHAGEL LARGE V-CAR-APS
502-812	AR-2092	EDTA, 50 G BOTTLE
502-820	AR-568	NICKEL CAR-APSULE 7.0 x 12 mm FB
502-822	AR-568	NICKEL CAR-APSULE 7.0 x 12 mm FB
502-825	AR-825	TIN CAR-APSULES 9MM X 21MM
502-834	AR-0	PLEASE CONTACT ALPHA FOR ENTIR
502-835	AR-0	PLEASE CONTACT ALPHA FOR ENTIR
502-844	AR-2304	COPPER STICKS, RESEALABLE BOTT
502-878	AR-189	COPPER STICKS 100 G AMPOULE
503-007	AR-007	DETACHABLE FLOAT VALVE
503-020	AR-3020	DISCONTINUED
503-031	AR-905	THREE-LEGGED FUNNEL
503-032	AR-1618	TWEEZERS, STAINLESS STEEL, CUR
503-032	AR-901	Regular Glass Scoop .14cc Vol.
503-033	AR-033	SULFUR TRAR-AP, FILLED
503-051	AR-051	OXYGEN BAFFLE 1 1/8" TUBE
503-052	AR-052	BAFFLE 1 1/4" TUBE
507-010	AR-911H	CATALYST TUBE FILLED WITH HOPC
507-011	AR-911E	CATALYST TUBE, UNFILLED
507-012	AR-012	REDUCER
507-100	AR-7100	CATALYST HEATER REPAIR
507-200	AR-7200	CATALYST HEATER REPAIR
507-207	AR-911	CATALYST TUBE, FILLED WITH PLA
507-207	AR-911R	CATALYST TUBE FILLED WITH RARE
507-300	AR-7300	CATALYST HEATER REPAIR
507-400	AR-7400	CATALYST HEATER REPAIR
507-500	AR-7500	CATALYST HEATER REPAIR
507-600	AR-7600	CATALYST HEATER REPAIR
507-700	AR-7700	CATALYST HEATER REPAIR
507-800	AR-7800	CATALYST HEATER REPAIR
507-900	AR-7900	CATALYST HEATER REPAIR
508-300	AR-8300	CATALYST HEATER REPAIR

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
508-400	AR-8400	CATALYST HEATER REPAIR
508-500	AR-8500	CATALYST HEATER REPAIR
508-600	AR-8600	CATALYST HEATER REPAIR, WITHOU
511-036	AR-182	EVAC PIN TUBE 150MM X 4MM ID,
511-038	AR-1038	Pin Tube Holder 3.5mm & 4.0mm
511-039	AR-181	EVAC PIN TUBE 150MM X 3.5MM ID
511-390	AR-183	EVAC PIN TUBE 150MM X 6MM ID,
511-390	AR-183M	REPLACED BY AR183
512-651	AR-2651	ABRAS. DISC 13 3/4 1 3/4
516-000	AR-906	PURIFYING TRAIN
516-004	AR-504	ROTOMETER
516-005	AR-505	DRYING TOWER
516-006	AR-506	ACID TOWER FOR H2SO4
516-063	AR-507	ITEM DISCONTINUED - NO REPLACE
517-204	AR-204	RUBBER SQUEEZE BULB
518-016	AR-016	POTASSIUM IODATE BURETTE
518-024	AR-024	LIGHT DIRECTOR ~;~
518-030	AR-017	DISCONTINUED - NO REPLACEMENT
518-045	AR-045	DISCONTINUED
519-004	AR-004	HYDROCARBON SULFUR TUBE
519-005	AR-900	IGNITOR
519-009	AR-009	PEDESTAL (2.25 INCHES)
521-084	AR-1084	VARIAC ~;~
521-112	AR-114	SILICONE TUBING 3/8" ID, 4 INC
521-113	AR-113	SILICONE RING SEAL
521-186	AR-187	SOLID STATE RECTIFIER 500 pf
521-254	AR-254	POWER RELAY (R) ~;~
521-287	AR-287	SUPPRESSOR COIL ~;~
521-305	AR-932	DRYING TUBE, 5.0"
528-018	AR-3818	CERAMIC CRUCIBLE, BAGGED, BOX
528-018	AR-3818F	CERAMIC CRUCIBLE, BOX OF 1000,
528-018	AR-8028F	FILTERING CRUCIBLE, SMALL, FOI
528-018-500	AR-3818F-500	CERAMIC CRUCIBLE, BOX OF 500,
528-028	AR-8028	FILTERING CRUCIBLE, SMALL
528-028	AR-8028F	FILTERING CRUCIBLE, SMALL, FOI
528-030	AR-8030	FILTERING CRUCIBLE, LARGE, 100
528-040	AR-040	COVER 4MM HOLE 1000/BOX
528-042	AR-042	POROUS COVER 1000/BOX
528-043	AR-9880	CER. COVER, 10 MM HOLE, 1000/BO
528-043	AR-9880S	10 MM HOLE COVER, BAG OF 250
528-050	AR-3818B	CERAMIC CRUCIBLES, BOX OF 500
528-053	AR-8053	COMBUSTION BOAT, BOX OF 500
528-053	AR-8053D	COMBUSTION BOAT, EACH
528-120	AR-8120	CERAMIC CRUCIBLE - LOW C & S,
528-203	AR-8203	COMBUSTION BOAT, BOX OF 500
528-203-250	AR-8203-250	COMBUSTION BOAT, BOX OF 250
528-230	AR-5282	GLAZED COMBUSTION BOATS

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
528-313	AR-8313	CERAMIC TUBE OBE 1 1/8 ID X 30
528-314	AR-8314	CERAMIC TUBE OBE 1 1/4 ID X 3
528-331	AR-8331	CERAMIC TUBE TE 1 1/8 ID X 27
528-332	AR-8332	CERAMIC TUBE TE 1 1/8 ID X 30
528-333	AR-8333	CERAMIC TUBE TE 1 1/4 ID X 30
529-007	AR-529	CERAMIC PEDESTAL, 3.05" "
529-011	AR-9011	CERAMIC LINER SC32/132
529-012	AR-9012	BOAT STOP KIT, 2 BOAT STOPS
529-018	AR-3203	ZIRCONIUM NOZZLE
529-022	AR-9022	BOAT STOP HALF SC32/132
529-026	AR-026	CRUCIBLE FOR FP 428 AND CHN-60
529-028	AR-9023	BOAT STOP, COARSE
529-031	AR-9031	CONE SUPPORT D
529-042	AR-9042	CERAMIC CRUCIBLE GLAZED TGA, 1
529-043	AR-9043	CERAMIC COVER TGA, FOR USE WIT
529-046	AR-9046	CERAMIC TUBE TF-10
529-047	AR-9047	CERAMIC CRUCIBLE GLAZED TGA,20
529-048	AR-9048	CERAMIC COVER TGA, FOR USE WIT
529-054	AR-9023	BOAT STOP, COARSE
529-055	AR-026	CRUCIBLE FOR FP 428 AND CHN-60
529-203	AR-8203	COMBUSTION BOAT, BOX OF 500
529-203	AR-8203P	COMB. BOATS SC32/132/CR12, BAG
529-204	AR-9204	COMBUSTION BOAT, S144, BOX OF
529-204	AR-9204-250	COMBUSTION BOAT, CS144, 250 PE
529-204	AR-9204P	COMBUSTION BOAT, S144, BAG OF
529-204-250	AR-9204-250	COMBUSTION BOAT, CS144, 250 PE
532-003	AR-535	HOUSING SEAT ASSEMBLY ~;~
532-060	AR-045	DISCONTINUED
532-065	AR-065	BURET FOR 532-500 ~;~
532-066	AR-066	VESSEL FOR 532-500 ~;~
532-300	AR-003	VALVE HOUSING AND SEAT
532-500	AR-532	AUTOMATIC SULFUR TITRATOR
533-104	AR-801	STARCH DISPENSING PIPETTE
537-103	AR-103	POLYFLOW TUBE 1/4" PER FT
540-140	AR-140	BULB ~;~
544-141	AR-141	THERMOMETER WELL
544-147	AR-4147	LEVELING BOTTLE
549-007	AR-907	ARMATURE ~;~
549-072	AR-9072	RESTRICTOR
549-076	AR-176	BULB, SOURCE FOR 518,532
550-004	AR-404	SUGGEST AR404S OR AR5210 REBUI
550-004	AR-404S	OSCILLATOR TUBE, SURPLUS
550-007	AR-5007	TIME DELAY TUBE DISCONTINUED
550-044	AR-044	DISCONTINUED - NO REPLACEMENT
550-046	AR-046	RESISTOR, 521/523
550-047	AR-0047	RF CHOKE/521/523 : DISCONTINU
550-053	AR-0053	POWER SWITCH FOR 521/523 ~;~

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
550-119	AR-019	RIGHT ANGLE BEND
550-120	AR-120	SULFUR COMBUSTION TUBE
550-122	AR-122	CARBON COMBUSTION TUBE
550-123	AR-123	INSERT, PEDESTAL
550-124	AR-124	PEDESTAL, 2.50"
550-162	AR-162	LOWER RING ADAR-APTER ~;~
550-165	AR-165	SPRING CLIP FOR COMB TUBE
550-167	AR-550	FURNACE FAN, 521/523
550-182	AR-180	QUARTZ ENCLOSED GRAR-APHITE ~;~
550-183	AR-179	QUARTZ CRUCIBLE INSERT
550-184	AR-184	DISCONTINUED, SUGGEST AR 263 A
550-187	AR-186	CAR-APACITOR - SURPLUS ITEM
563-035	AR-035	DISCONTINUED REPL BY NOTHING ~;~
563-037	AR-037	DISCONTINUED
566-008	AR-006	INDUCTION COIL
568-016	AR-428	O RING
572-016	AR-72016	CONTROL VALVE FOR 572
572-018	AR-572	DISCONTINUED ~;~
572-029	AR-029	FLOAT VALVE
572-031	AR-2031	SPLASH PROOF DEVICE
578-672	AR-672	"O" RING
585-235	AR-200	PIN CUTTER WHEEL , 7" DIAMETER
585-270	AR-270	PIN CUTTER WHEEL , Al ₂ O ₃ 7 I
585-359	AR-5359	SPACER FOR AR2000, 1/2 GRAM
585-360	AR-5360	SPACER FOR AR2000, 1 GRAM
585-600	AR-2000	*PIN CUTTER, WITH SPACERS AND
589-305	AR-9305	1/4" RUBBER SLEEVE
589-364	AR-9364	OVEN FAN ~;~
589-551	AR-422	O RING (AUTO DROP)
589-699	AR-699	LANCE SLEEVE, RUBBER
589-859	AR-859	NOT CURRENTLY AVAILABLE ~;~
589-887	AR-887	O - RING
589-891	AR-9891	TIMER RELAY 750/589 ~;~
593-100	AR-593	TIMER, SULFUR TITRATOR
598-142	AR-859	NOT CURRENTLY AVAILABLE ~;~
598-157	AR-157	H ₂ O TRAR-AP FILLED 750,589
598-158	AR-158	H ₂ O TRAR-AP UNFILLED
598-257	AR-257	DISCONTINUED
598-614	AR-613	BULB/SWITCH 750-100 6.3V ~;~
601-064	AR-064	QUARTZ HEATER TUBE, CF-10
601-065	AR-1065	PYREX REAGENT TUBE, CF-10
601-132	AR-132	COMB TUBE FOR CNH600/800
601-133	AR-1133	"O" RING
601-158	AR-1158	O RING, CHN 900
601-159	AR-1159	O RING, CHN 900
601-201	AR-1201	REACTION TUBE
601-390	AR-390	QUICK DISC TUBE

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
601-397	AR-397	CERAMIC WOOL SHIM, PACKAGE OF
601-437	AR-437	QUARTZ CATALYST TUBE FP428
601-438	AR-438	COMBUSTION TUBE, FP428, INCLUD
601-441	AR-441	QUARTZ LANCE TUBE, FP428
601-442	AR-442	O RING, FP428
601-482	AR-482	HONEYCOMB SUPPORT, FP428
601-483	AR-483	PRINTER PAR-APER 9 1/2 x 11" , 600
601-484	AR-484	RIBBON SPOOL PRINTER
601-490	AR-1490	FOILER FIXTURE
601-504	AR-1504	O RING, LOADING BLOCK FP228
601-505	AR-1505	"O" RING, FP/NS 2000
601-506	AR-1506	COMBUSTION TUBE, FP228
601-508	AR-438	COMBUSTION TUBE, FP428, INCLUD
601-654	AR-1654	THERMOCOUPLE FOR CHN2000/CHN10
601-684	AR-684	O RING
601-691	AR-1691	.125 CFLEX TUBING PER FT
601-691-030	AR-1691-30	CFLEX TUBE, 30" SECTION
601-691-036	AR-1691-36	CFLEX TUBING, .125" ID
601-862-110	AR-1862	DISCONTINUED, NO REPLACEMENT
601-920	AR-920	O RING
601-939	AR-1939	O-RING DOSING VALVE
601-963	AR-AED4032	SMOOTH WALL TIN CAR-APSULES FLAT
601-963	AR-1963	DISCONTINUED
601-981	AR-981	SIDE ARM CLEANING BRUSH
602-197	AR-6186	FTG ADAR-APTER, RC612
602-615	AR-2615	CERAMIC FURNACE HOUSING ~;~
602-802	AR-802	O RING, FOR RH 402
602-803	AR-803	O RING, FOR RH 402
602-969	AR-969	5000 DETERMINATIONS PKG ~;~
602-977	AR-AEE2374	SIDE ARM PLUG
603-047	AR-63047	URETH TUBING, .061 X ~;~
603-048	AR-63048	URETH TUBING, .061 NAVY BLUE
603-049	AR-63049	URETH TUBING, .061 X ~;~
603-050	AR-63050	URETH TUBING, BLACK .061
603-051	AR-63051	URETH TUBING, WHITE .061
603-052	AR-63052	URETH TUBING, .061 X .12 ~;~
603-101-460	AR-980	IN LINE PARTICLE FILTER
603-350	AR-AEE2464	O RING INSTALLATION TUBE
603-433	AR-3433	O-RING
603-493	AR-3493	O RING
603-706	AR-706	PEDESTAL RH402
603-707	AR-707	REACTION TUBE RH 402 ~;~
603-710	AR-710	QUARTZ CRUCIBLE PACK/24
603-910	AR-3910	O RING
603-970	AR-3970	O-RING, LOADING MECHANISM, RH4
603-971	AR-3971	O-RING, LOADING MECHANISM FOR
604-112	AR-604	QUARTZ CATALYST TUBE ~;~

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
604-114	AR-923	PEDESTAL ONLY
604-163	AR-163	BRASS BRUSH
604-172	AR-172	O RING
604-241	AR-4241	REACTION TUBE, RH 402
604-271	AR-271	SS BRUSH DIAGONAL SLOT
604-370	AR-AEE2930	FILAMENT MATCHED PAIR
604-373	AR-4373	CAR-APSULE HOLDER, FOR USE WITH A
604-496	AR-4496	TEFLON TUBING 1/8" O.D., FOOT
604-518	AR-518	SS BRUSH, NO SLOT, FINE WIRE
604-541	AR-4541	LANCE TUBE, CHN900/932, QUARTZ
605-520-087	AR-6055	O RING PACKAGE
605-814	AR-5814	PEDESTAL ADAR-APTER ASSY, INCLUDE
605-815	AR-5815	PEDESTAL ADAR-APTER, FOR USE WITH
605-816	AR-816	CERAMIC PEDESTAL
605-817	AR-817	PEDESTAL ADAR-APTER, CS125 ~;~
605-818	AR-818	PEDESTAL ADAR-APTER ~;~
605-828	AR-828	O RING, FP428
605-830	AR-830	O-RING
605-875	AR-5875	FURNACE RED/OX TUBE
606-026	AR-6026	GRAR-APHITE CRUCIBLE, PKG 10
606-202	AR-202	FURNACE FAN FP2000 SERIES
606-224	AR-6224	THERMOCOUPLE
606-286	AR-563	SHORT TERMINAL BRAIDED STRAR-AP,
606-308	AR-6308	STEPPED CERAMIC BOAT STOP
606-309	AR-6309	OUTER COMBUSTION TUBE
606-310	AR-1477	O RING, OUTER COMBUSTION TUBE
606-310	AR-6310	REPLACED BY AR1477 ~;~
606-312	AR-1476	O RING, INNER COMBUSTION TUBE
606-312	AR-6312	REPLACED BY AR1476 ~;~
606-313	AR-6313	INNER COMBUSTION TUBE
606-316	AR-6316	ADAR-APTOR ASSEMBLY, SC 432, WITH
606-317	AR-6317	PURGE TUBE, WITH FITTINGS
606-318	AR-6318	3 HEATING ELEMENTS SC432
606-319	AR-6319	BRAIDED LEAD WIRE, 6.5 INCHES
606-320	AR-6320	BRAID CLAMP FOR ELEMENT
606-327	AR-6327	GAS DRYING TUBE, SC-432
606-333	AR-6333	O-RING, LANCE SC32
606-340	AR-6340	LANCE TUBE FP2000 250MM
606-403	AR-403	COLLECTION TUBE BRUSH
606-444	AR-6644	# 2 STEPPED CERAMIC STOP
606-547	AR-6547	BOAT PULLER, SC 432
606-550	AR-6550	FINGER ARM FOR AR6547 BOAT PUL
606-601	AR-6601	HEATING ELEMENTS
606-601-100	AR-6601-100	HEATING ELEMENT KIT
606-602	AR-6602	FURNACE END CAR-APS, FP2000
606-655	AR-6655	LARGE VOLUME DRYING TUBE ~;~
606-747	AR-6747	DRILLED LANCE TUBE, 2HOLE ~;~

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
607-073	AR-607	O RING
608-065	AR-8065	LANCE TUBE TRUSPEC MICOR 02
608-089	AR-8089	O-RING
608-091	AR-8091	O-RING
608-130	AR-AEC7000	INNER PYROLYSIS TUBE
608-339	AR-8339	CHN-1000 LANCE TUBE
608-342	AR-342	INNER COMB TUBE, QUARTZ
608-343	AR-8343	HEATING ELEMENT, CHN-1000
608-344	AR-344	OUTER COMB TUBE, CERAMIC
608-378	AR-1476	O RING, INNER COMBUSTION TUBE
608-379	AR-8379	QUARTZ WOOL STRIPS, 15"
608-505	AR-8505	INNER COMB TUBE, CERAMIC
608-617	AR-AP149C	IN-LINE FILTER, FOR PARR BOMB
608-624	AR-624	O-Ring
608-866	AR-866	.1875 ID CFLEX TUBING, FT
609-107	AR-241	SILICONE GREASE 150g
609-114	AR-9114	COMB. TUBE O-RING FP428
609-116	AR-9116	COMB TUBE NEW STYLE, FP 428, W
609-129	AR-9129	LANCE, STRAIGHT
609-170	AR-9170	QUICK CAR-AP KIT INCLUDES ~;~
609-171	AR-9171	PLUNGER ASSEMBLY ~;~
609-172	AR-9172	TUBING, CELLULOSIC, .210 ID, 2
609-271	AR-9271	RADIAL SHAFT SEAL, FOR FP2000
609-394	AR-394	HEAT SHIELD, FP/NS2000
609-474	AR-9474	FINGER ARM FOR AR1207 PULLER A
609-556	AR-AEE2194	MOTOR
609-728	AR-9728	METAL SHAFT PUSH/PULL BEARING
610-735	AR-735	SS BRUSH, DIAGONAL SLOT, COARS
611-207	AR-1207	BOAT PULLER ASSEMBLY, FP2000,
611-255	AR-9171	PLUNGER ASSEMBLY ~;~
611-351-142	AR-1142	CATALYST TUBE (RHEN602)
611-351-142	AR-279	CATALYST HEATER TUBE 10.75"
611-351-142	AR-9412	CATALYST TUBE, 10.75 INCH, TC-
611-351-180	AR-51180	ELECTRODE INSERT RHEN602 WITH
611-351-181	AR-691	INSERT ELECTRODE TC/OHN ANALYZ
611-351-182	AR-691	INSERT ELECTRODE TC/OHN ANALYZ
611-351-183	AR-51183	ELECTRODE INSERT RHEN602 WITHO
611-351-223	AR-1351	COMBUSTION TUBE, QUARTZ, WC600
611-431	AR-1431	O-RING FOR FP2000 LANCE (ORANG
611-476	AR-1476	O RING, INNER COMBUSTION TUBE
611-477	AR-1477	O RING, OUTER COMBUSTION TUBE
611-790	AR-790	"O" RING
612-340	AR-2340	O RING
613-563	AR-563	SHORT TERMINAL BRAIDED STRAR-AP,
614-591	AR-257	DISCONTINUED
614-691-110	AR-4961	LARGE CRUCIBLE FOR SOLID SAMPL
614-822-101	AR-4101	SMALL BOAT LOOP, NICKEL

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
614-822-102	AR-4181	Small Nickel Boat
614-822-102	AR-4822	Mercury Catalyst Tube
614-822-113	AR-4113	LARGE BOAT LOOP, NICKEL
614-822-114	AR-4114	Large Nickel Boat
614-961	AR-4961	LARGE CRUCIBLE FOR SOLID SAMPL
614-961-110	AR-4961	LARGE CRUCIBLE FOR SOLID SAMPL
614-962	AR-4962	LANCE TUBE, FP 528
614-963	AR-AEE2498	LANCE ASSEMBLY COMPLETE
616-063	AR-6063	BOAT PULLER ASSEMBLY
616-068	AR-6068	O-RING, INNER COMB TUBE, S-144
616-070	AR-6070	INNER COMBUSTION TUBE
616-086	AR-6086	STEPPED BOAT STOP
616-089	AR-6089	OUTER COMBUSTION TUBE
616-108	AR-6108	COMBUSTION TUBE FOR FP528
616-138	AR-6138	O-RING
616-139	AR-6139	O-RING, FP-528
616-146	AR-6146	LANCE, FP 528
616-152	AR-6152	CRUCIBLE EXTRACTOR TOOL, FP 52
616-504	AR-6504	PARTICLE FILTER
616-540	AR-AEE2499	DOSING VALVE STEM
616-622	AR-6622	MICRO SWITCH CN 2000 SWITCH ON
616-632	AR-6632	LANCE TUBE 253MM, FP2000 & SC
616-635	AR-394	HEAT SHIELD, FP/NS2000
617-031	AR-7031	O-RING
617-083	AR-7083	O-RING
617-136	AR-7136	O RING
617-154	AR-3930	ELECTRODE INSERT REPAIR, CUSTO
617-440	AR-440	O-RING
617-441	AR-7441	MICRO SCREEN, TCH600
617-548	AR-7548	REAGENT TUBE, 11.4 INCH
617-605	AR-7605	LARGE CRUCIBLE FOR LIQUID, FP
617-716	AR-7716	REAGENT TUBE
618-376	AR-376	ELECTRODE INSERT TIP USED WITH
618-482	AR-691	INSERT ELECTRODE TC/OHN ANALYZ
618-603	AR-603	CATALYST TUBE, 13.0 INCH
618-657	AR-6186	FTG ADAR-APTER, RC612
618-809	AR-8809	REAGENT TUBE, 11.4 INCHES
619-065	AR-9065	COMBUSTION TUBE FOR TRUSPEC
619-154	AR-9154	CATALYST HEATER TUBE
619-192	AR-9192	O RING TruSpec Loading Head
619-233	AR-9233	THERMOCOUPLE 90 DEGREE BEND
619-268	AR-9268	HIGH SENSITIVITY REAGENT TUBE
619-316	AR-9316	REAGENT TUBE
619-425	AR-9425	STRAIGHT GLASS TUBE, TRUSPEC
619-455	AR-9455	LARGE REAGENT TUBE
619-566	AR-9566	BRAIDED LEAD WIRE
619-590-409	AR-0409	QUARTZ REAGENT TUBE CS744 CS84

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
619-590-412	AR-0412	QUARTZ CATALYST TUBE CS744 CS8
619-590-861	AR-861	O RING CS744 CS844
619-591-314	AR-1314	REAGENT TUBE CS744/S844 OHN836
619-591-315	AR-1315	FURNACE REAGENT TUBE CS744 CS8
619-591-385	AR-1385	O RING
619-591-683	AR-6195	NO LONGER AVAILABLE 03/05/2014
619-591-683	AR-683	COMBUSTION TUBE FOR CS744 AND
619-591-699	AR-1699	PARTICLE FILTER
619-591-740	AR-1740	REAGENT TUBE CS744 CS844
619-591-755	AR-920	O RING
619-591-969	AR-1969	BRASS BRUSH
619-591-970	AR-1970	SLOTTED SS CLEANING BRUSH
619-592-403	AR-291	OXYGEN REGULATOR
619-592-926	AR-920	O RING
619-774	AR-9774	LANCE, TRUSPEC, CHNS MICRO
619-776	AR-9776	REAGENT TUBE CHNS 233 MM
619-777	AR-9777	DIRECTING MIDDLE TUBE CHNS
619-803	AR-3302	BRUSH END MODIFIED
619-880	AR-9880	CER. COVER, 10 MM HOLE, 1000/BO
619-880	AR-9880S	10 MM HOLE COVER, BAG OF 250
619-880-100	AR-9880-100	CERAMIC COVER, 10MM HOLE, 100/
619-895	AR-619	GRAR-APHITE CRUCIBLE, RHN600, BOX
619-895	AR-6198	ELECTRODE INSERT FOR 619-895 (
619-895	AR-9895	GRAR-APHITE DOSING CRUCIBLE, BOX
619-895-100	AR-619-100	GRAR-APHITE CRUCIBLE, RHN600
619-896	AR-6198	ELECTRODE INSERT FOR 619-895 (
619-926	AR-3930	ELECTRODE INSERT REPAIR, CUSTO
619-926	AR-692	LOWER ELECTRODE MALE PLATE
620-814	AR-691	INSERT ELECTRODE TC/OHN ANALYZ
620-976	AR-3930	ELECTRODE INSERT REPAIR, CUSTO
620-985	AR-3930	ELECTRODE INSERT REPAIR, CUSTO
621-331	AR-9047	CERAMIC CRUCIBLE GLAZED TGA,20
622-001-026	AR-1026	GLASS TUBE FOR 628 SERIES
625-067	AR-5067	HEATING ELEMENT - 6 RODS
625-068	AR-25068	CLAMP FOR AR5067 ELEMENTS 12PK
625-235	AR-816	CERAMIC PEDESTAL
625-401-294	AR-1317	STAINLESS STEEL PURGE TUBE
625-401-301	AR-01301	ORING
625-401-302	AR-01302	O RING
625-401-302	AR-1302	ORING SC632/832
625-401-317	AR-1317	STAINLESS STEEL PURGE TUBE
625-401-351	AR-6254	LANCE TUBE SC632 W/AUTLOADER
625-401-352	AR-1352	LANCE TUBE, TRUSPEC, 244 MM
625-501	AR-6255	THERMOCOUPLE, RC612
625-501-358	AR-1358	THERMOCOUPLE FOR SC632
625-505-125	AR-5125	CATALYST HEATER TUBE, RC612
625-505-138	AR-5138	FURNACE HEATING ELEMENT RC612

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
625-505-168	AR-5168	COMBUSTION TUBE, RC612
625-505-271	AR-5271	KNURLED ADAR-APTER, RC612
625-510-402	AR-5104	COMBUSTION TUBE, RC 612
625-510-566	AR-3302H	HANDLE FOR AR3302 ELECTRODE BR
625-510-614	AR-5271	KNURLED ADAR-APTER, RC612
625-511-154	AR-511154	LOWER CRUCIBLE, PKG OF 5
625-601-554	AR-1554	OUTER COMBUSTION TUBE
625-601-555	AR-1555	INNER COMBUSTION TUBE
625-601-904	AR-01904	O-RING
625-602-187	AR-602	CERAMIC LANCE TUBE USED W/AUTO
625-602-297	AR-02297	O-RING
625-602-539	AR-6327	GAS DRYING TUBE, SC-432
625-710-392	AR-0392	REAGENT TUBE OHN836
625-711-116	AR-1116	REAGENT TUBE OHN836
625-711-117	AR-1117	REAGENT TUBE OHN836
625-711-182	AR-11182	UPPER ELECTRODE CLEANING BRUSH
625-711-186	AR-11186	UPPER ELECTRODE CLEANING BRUSH
625-711-599	AR-5933	COOLING AGENT, 1 PINT CONTAINE
626-001-691	AR-01691	CERAMIC LANCE NO AUTLOADER
626-001-692	AR-01692	PURGE TUBE W/O AUTLOADER
626-001-871	AR-01871	GAS DRYING TUBE SC832
626-002-467	AR-02467	PURGE TUBE USED W/AUTOLOADER
633-101-361	AR-563	SHORT TERMINAL BRAIDED STRAR-AP,
633-101-362	AR-9566	BRAIDED LEAD WIRE
633-101-363	AR-9566	BRAIDED LEAD WIRE
633-101-460	AR-980	IN LINE PARTICLE FILTER
633-103-225	AR-7548	REAGENT TUBE, 11.4 INCH
660-010-130	AR-9233	THERMOCOUPLE 90 DEGREE BEND
701-030	AR-APT1030	TORCH FOR LECO ICP2000 PN: 701
702-557	AR-2557	O RING
704-295-109	AR-1585	OVER TEMP THERMOCUPLE, RC612
704-295-109	AR-95109	THERMOCOUPLE
722-522	AR-7073	THERMOCOUPLE FOR USE ON SC32/1
722-568	AR-233Q	CERAMIC ADHESIVE
736-026	AR-027B	ALPHACEL III -20+40 MESH, 30 P
750-301	AR-933	DRYING TUBE, 9.0"
751-300	AR-1300	NO LONGER VALID - PLEASE SEE A
760-034	AR-034	GRAR-APHITE CRUCIBLE, BOX OF 1000
760-072	AR-072	28 VDC BULB
760-138	AR-138	TWEEZER, STRAIGHT, SELF-CLOSIN
760-139	AR-766	DELUXE CRUCIBLE TWEEZERS, FOR
760-144	AR-022B	CLEANING BRUSH /90 BEND
760-145	AR-0145	FUNNEL FOR POWDER SAMPLE, USE
760-193	AR-193	METAL SCREEN GLASS TUBE
760-224	AR-224	O RING
760-225	AR-425	"O" RING (HEAD PLUNGER)
760-226	AR-428	O RING

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
760-293	AR-293	BRUSH, CLEANING, SMALL
760-301	AR-931	DRYING TUBE, 1.5"
760-304	AR-0304	TEFLON WASHER
760-316	AR-316	WASHER, FLAT - QUICK DISC ~;~
760-317	AR-317	GLASS TUBE 6.64"
760-333	AR-333	GLASS TUBE, 8.0" OR 203 MM
760-335	AR-424	"O" RING (HEAD TOP PLATE)
760-359	AR-359	PLUNGER STEM O-RING
760-414	AR-414	GRAR-APHITE CRUCIBLE, BOX OF 1000
760-414	AR-414P	GRAR-APHITE CRUCIBLE
760-841	AR-194	"O" RING QUICK DISCONNECT
761-047	AR-047	3" DRYING TUBE
761-133	AR-133	DRYING TUBE, 17.125"
761-190	AR-190	GLASS TUBE 11.75" OR 298 MM
761-383	AR-293H	BRUSH HANDLE, KNURLED
761-492	AR-492	"O" RING, FP428
761-618	AR-618	DISCONTINUED
761-651	AR-1609	"O" RING, LANCE, FP/CHN
761-739	AR-639	TIN FLUX 0.5GM PELLETS, 1 POUN
761-747	AR-611	SCHUTZE REAGENT, 100 G PER BOT
761-814	AR-814	SUPPORT RING
761-918	AR-918	"O" - RING
761-929	AR-929	CRUCIBLE TONGS FOR AR3818
762-023	AR-023	NUT ADAR-APTER SEAL
762-056	AR-056	COMBUSTION TUBE
762-057	AR-416	"O" RING (DUST FILTER)
762-058	AR-058	O RING (LOWER ADAR-APTER)
762-059	AR-614	"O" RING (CUP ASSEMBLY)
762-060	AR-417	O RING
762-062	AR-913	OXYGEN LANCE
762-066	AR-068	RECTIFIER (R) ~;~
762-066	AR-7620	NO LONGER AVAILABLE!!
762-192	AR-192	DISCONTINUED
762-357	AR-357	"O" RING (UPPER ADAR-APTER)
762-501	AR-2501	HANDLE FOR REBUILT 521 ~;~
762-621	AR-2621	IR CO2 DETECTOR CS44,IR12
762-627	AR-2627	BISCUIT FAN FOR CS44
762-695	AR-695	Tin Metal Accelerator 20 Mesh
762-697	AR-697	BRUSH (NO HANDLE)
762-741	AR-650	HON IN TITANIUM STANDARD H=0.0
762-747	AR-546	HYDROGEN IN STEEL PIN STANDARD
762-823	AR-697H	KNURLED HANDLE/AR 697
762-936	AR-936	V2O5 Glass Scoop 0.6cc Vol
763-027	AR-266B	ALPHACEL TUNGSTEN BULK, 30 POU
763-028-500	AR-027V	ALPHACEL III TUNGSTEN, 20-40 M
763-029	AR-2345	NICKEL BASKETS 1.5G BASKETS, 1
763-065	AR-2344	NI BASKETS LOW OXYGEN, 1g, 10

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
763-191	AR-103	POLYFLOW TUBE 1/4" PER FT
763-213	AR-213	GRAR-APHITE CRUCIBLE RH-1E, 1000
763-213	AR-213P	GRAR-APHITE CRUCIBLE RH-1E, 100 P
763-261	AR-261	FILTER PAR-APER .510" " OD, PACK O
763-263	AR-266	ALPHACEL TUNGSTEN -12+20
763-265	AR-265	GLASS WOOL HALF POUND
763-266	AR-027	ALPHACEL TUNGSTEN III -20+40 M
763-266	AR-7629	ALPHACEL V, -20/+40 MESH TUNGS
763-295	AR-295	THERMOSTAT/UPPR ELECTRODE ~;~
763-379	AR-3379	O RING, FP428
763-467	AR-467	IRON CHIP , 25 POUNDS BULK OR
763-467	AR-467C	DISCONTINUED
763-649	AR-749	RED LED LAMP MV5023 ~;~
763-873	AR-3873	BULB, ISOLATION BOX
763-930	AR-930	ACTIVATED CHARCOAL - ITEM DISC
764-015	AR-015	BULB, 532-500, 28V VDC ~;~
764-078	AR-4078	DISCONTINUED
764-098	AR-098	CELLULOSE TRAR-AP MATERIAL 1lb
764-215	AR-215	REPLACED BY AR4216
764-216	AR-4216	INERT GAS REGULATOR
764-242	AR-242	ALPHACEL TIN PELLETS, 1 POUND
764-330	AR-331	GRAR-APHITE CRUCIBLE - BAG of 100
764-330	AR-331S	GRAR-APHITE CRUCIBLE .625 OD, SOL
764-490	AR-490	PEDESTAL 2-3/4" " HIGH ~;~
764-502	AR-4544	ROTRON FURNACE FAN ~;~
764-545	AR-2832	LUBE OIL, .25%S , 100ml
764-546	AR-2835	LUBE OIL, 1.00%S , 100ml
764-547	AR-2836	LUBE OIL, 2.00%S, 100 ML BOTTL
764-640	AR-334	100 PF CAR-APACITOR
764-641	AR-335	200 PF CAR-APACITOR
764-720	AR-064	QUARTZ HEATER TUBE, CF-10
764-762	AR-762	GLASS TUBE 11.88"
764-775	AR-775	WR-12 TRAR-AP MATERIAL
765-068	AR-5068	DISCONTINUED
765-523	AR-523	OVEN FAN WITH ADAR-APTER
765-717	AR-917	OXYGEN LANCE (WR-12)
765-746	AR-746	ROTOMETER, WR12 ~;~
765-818	AR-5818	OPTIC CABLE ~;~
765-847	AR-847	HEAT EXCHANGE FAN TC 136 ~;~
765-933	AR-5933	COOLING AGENT, 1 PINT CONTAIN
765-976	AR-5976	O RING FOR HEATED TUBES BVN
766-053	AR-766	DELUXE CRUCIBLE TWEEZERS, FOR
766-363	AR-363	DISCONTINUED NO REPLACEMENT
766-363	AR-363P	GRAR-APHITE CRUCIBLE RH-1 ~;~
766-557	AR-6557	O RING ~;~
766-646	AR-6646	COATING AIR FILTER 10 OZ
766-823	AR-823	MALE CLIPLOCK FITTING ~;~

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
767-243	AR-243	FEMALE CLIPLOCK FITTING
767-277	AR-277	GRAR-APHITE CRUCIBLE, BOX OF 1000
767-277	AR-277P	GRAR-APHITE CRUCIBLE, BAG OF 200
767-385	AR-085	QUARTZ WOOL 12+ μ , COARSE WADD
767-473	AR-473	UTILITY FUNNEL FOR GLASS TUBES
767-500	AR-914	OXYGEN LANCE (IR 33)
767-518	AR-7518	RETAINING SPRING CLIP
767-541	AR-541	DRYING TUBE, BUBBLE, 7 INCHES
767-542	AR-542	Tungsten Glass Scoop .1cc Vol
767-700	AR-9822	SARTORIUS BA 110 BALANCE ~;~
767-757	AR-1706	2.07% SULFUR COAL STANDARD
767-758	AR-1707	2.63% SULFUR IN COAL STANDARD
767-759	AR-720	COKE STANDARD, 0.77% SULFUR
767-759	AR-724	COKE STANDARD, 1.21% SULFUR
767-838	AR-838	""O"" - RING
768-774	AR-043S	REPLACED WITH AR9880 OR AR9880
768-785	AR-785	DISCONTINUED ~;~
768-798	AR-8798	CLEAN OUT PLUG
768-944	AR-944	""O"" RING
768-980	AR-980	IN LINE PARTICLE FILTER
769-310	AR-9310	O-RING
769-335	AR-469	RELAY
769-520	AR-9520	GRAR-APHITE CRUCIBLE, RH-2, PACKA
769-520	AR-9520C	GRAR-APHITE CRUCIBLE, RH-2, PKG O
769-520	AR-9520M	GRAR-APHITE CRUCIBLE, RH-2, P KG
769-569	AR-054	GRAR-APHITE CRUCIBLE LIDS, WITH H
769-608	AR-608	HALOGEN TRAR-AP Sb METAL
769-610	AR-610	HALOGEN TRAR-AP MATERIAL
769-761	AR-761	TALL GRAR-APHITE CRUCIBLE, ALUMI
769-821	AR-9821	DISCONTINUED
769-897	AR-897	DISCONTINUED UNTIL FURTHER NOT
769-898	AR-898	LOWER PED, RH-2 & RH-3
769-994	AR-414	GRAR-APHITE CRUCIBLE, BOX OF 1000
769-994	AR-994	REPLACED BY AR 414
770-247	AR-247	FAN FILTER, PACK OF 5
770-436	AR-523	OVEN FAN WITH ADAR-APTER
770-718	AR-0718	DISCONTINUED!!! PRINTER
770-764	AR-764	PRINTER PAR-APER ROLL
770-778	AR-778	RIBBON SPOOL ASSEMBLY
770-779	AR-7707	NO LONGER AVAILABLE!!!
771-015	AR-1015	POTENTIOMETER, CS46 ~;~
771-030	AR-030	BUBBLE TUBE
771-040	AR-3913	O RING
771-275	AR-1275	PNEUMATIC CYLINDER
771-354	AR-354	OSCILLATOR TUBE CS144-600, ELT
771-362	AR-1362	O-RING
771-367	AR-1367	QUARTZ CRUCIBLES

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
771-367	AR-1367S	QUARTZ CRUCIBLES
771-373	AR-373	DRYING TUBE, WITH SAF-T-WRAR-AP
771-401	AR-401	REACTION TUBE/RH-2, RH-3
771-407	AR-407	O RING
771-425	AR-1425	DISCONTINUED, NOT ABLE TO OFFE
771-446	AR-429	""O"" RING
771-606	AR-606	PEDESTAL ~;~
771-883	AR-083	BUBBLE TUBE, CATALYST
771-991	AR-991	O RING
772-005	AR-760	CHART PAR-APER ROLL, RH2 ~;~
772-269	AR-269	O RING
772-381	AR-2381	LANCE CLEANING WIRE, FOR CS DE
772-402	AR-2402	LANCE TUBE ASSEMBLY
772-403	AR-6403	CERAMIC LANCE TUBE 267mm
772-446	AR-446	VIEWING WINDOW
772-449	AR-449	VIEWING WINDOW ROD ~;~
772-456	AR-456	THIS IS NO LONGER AVAILABLE!!!
772-457	AR-456	THIS IS NO LONGER AVAILABLE!!!
772-460	AR-2460	PRESSURE TRANSDUCER
772-469	AR-469	RELAY
772-508	AR-962	HEATING ELEMENT
772-522	AR-7073	THERMOCOUPLE FOR USE ON SC32/1
772-560	AR-2560	CHOPPER MOTOR SC32
772-738	AR-738	O RING
772-739	AR-544	CATALYST TUBE, 10.75""
772-747	AR-2747	FAN
772-795	AR-795	OVEN FAN
772-905	AR-057	COMBUSTION TUBE,CS144-600
772-906	AR-2906	COMBUSTION TUBE SUPPORT
772-910	AR-2910	O RING,PISTON SEAL, FP228 (BLA
773-066	AR-3066	CERAMIC STANDOFF - Standoff; 8
773-076	AR-3076	DISCONTINUED, NO REPL.
773-106	AR-106	DELIVERY TUBE FOR SULFUR TITRA
773-136	AR-3136	FITTING, S.S.
773-140	AR-3140	NO LONGER AVAIALBLE!!!
773-141	AR-320	O RING
773-311	AR-311	SST FILTER
773-314	AR-3314	UNFLUTED CLEANING BRUSH
773-316	AR-3316	REPLACED BY AR789
773-316	AR-789	LANCE TUBE ASSEMBLY, 4.28 INCH
773-320	AR-320	O RING
773-392	AR-3932R	REPLACED BY AR3932
773-561	AR-561	SS BRUSH STR. SLOT, CS444
773-579	AR-579	Metal Scoop 0.18cc Volume
773-662	AR-3662	POLY FLO FLEX TUBING CHN
773-913	AR-3913	O RING
773-928	AR-3928	DELIVERY TUBE SC32/132

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
773-931	AR-9931	TEFLON DISCONNECT
773-932	AR-3932	COMBUSTION TUBE ASSEMBLY
773-935	AR-935	REPLACED BY AR 233
774-031	AR-4031	QUICK DISCONNECT TUBE,1.38OD
774-174	AR-4174	PURGE TUBE, FOR USE ON SC DETE
774-204	AR-4204	AC CRUCIBLE, STAINLESS STEEL,
774-205	AR-205	FIRING FUSE WIRE, 300 DETERMIN
774-208	AR-208	BENZOIC ACID, 1 GRAM, 10/PACK
774-208	AR-208C	BENZOIC ACID, 1 GRAM, 100/PAC
774-208	AR-208V	BENZOIC ACID TABLETS 1.0 GM,
774-211	AR-4211	"O"-RING 2.6X2.875X.005 ~;-
774-253	AR-4253	"O" RING
774-425	AR-432	O RING, CLOSURE BLOCK
774-460	AR-460	THERMOCOUPLE ASSEMBLY
774-489	AR-489	O RING, UPPER ELECTRODE (BROWN
774-570	AR-4570	VACUUM PUMP - 115V
774-571	AR-4571	VACUUM PUMP - 24VDC, SC-32CR-1
774-613	AR-4613	O-RING HF 400
774-638	AR-638	QUICK DISCO TUBE,6" STRT
774-649	AR-653	CLEANING BRUSH W/ HANDLE
774-651	AR-652	CLEANING BRUSH TC136
774-652	AR-652	CLEANING BRUSH TC136
774-690	AR-690	KNURLED RING COLLAR
774-691	AR-691	INSERT ELECTRODE TC/OHN ANALYZ
774-692	AR-692	LOWER ELECTRODE MALE PLATE
774-809	AR-809	INDUCTION COIL/HF10
774-841	AR-841	SENSOR, LOWER DRAIN
775-288	AR-3403	BENZOIC ACID POWDER, 454 G
775-306	AR-5306	10 MICRON FILTER, WITH AR520 O
775-312	AR-312	O-RING LOWER ELECTRODE
775-313	AR-313	TEFLON BACKUP RING ~;-
775-314	AR-314	O RING, PISTON BLOCK
775-330	AR-5330	FOUR WAY VALVE
775-355	AR-355	CATALYST HEATER TUBE 10.75 INC
775-416	AR-293	BRUSH, CLEANING, SMALL
775-416	AR-293H	BRUSH HANDLE, KNURLED
775-431	AR-431	REPLACED BY AR5431
775-431	AR-5431	INNER GRAR-APHITE CRUCIBLE, BOX O
775-431	AR-5431S	INNER GRAR-APHITE CRUCIBLE THICK
775-433	AR-433	OUTER CARBON CRUCIBLE 50 PER P
775-446	AR-3930	ELECTRODE INSERT REPAIR, CUSTO
775-446	AR-5446	DISCONTINUED - SEE AR3930 (NEE
775-601	AR-601	BUBBLE TUBE, LEFT, 6"
775-625	AR-625	PUSH BUTTON SWITCH
775-677	AR-677	O RING
775-722	AR-456	THIS IS NO LONGER AVAILABLE!!!
775-814	AR-402	REACTION TUBE/RH-2 THICK ~;-

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
775-859	AR-5859	WIRE CLEANING BRUSH
775-880	AR-880	GLASS DISCO TUBE,4" STRAIGHT.
775-892	AR-431P	INNER CARBON CRUCIBLE, 100 PER
776-036	AR-6036	PNEUMATIC (AIR) REGULATOR
776-095	AR-095	SILICONE TUBE, 19"
776-129	AR-129	INDUCTION COIL
776-234	AR-6234	10 MICRON FILTER, 1GPM
776-247	AR-053	REPLACED BY AR6247, DISCONTINU
776-247	AR-053HD	GRAR-APHITE CRUCIBLE, 1000/PACK
776-247	AR-053P	GRAR-APHITE CRUCIBLE, BAG OF 200
776-247	AR-6247	GRAR-APHITE CRUCIBLE, 1000/PACK
776-247	AR-6247M	GRAR-APHITE CRUCIBLE, RESEALABLE
776-247	AR-6247P	GRAR-APHITE CRUCIBLE, 200 pack
776-247	AR-625	PUSH BUTTON SWITCH
776-265	AR-6265	PEDESTAL POST ASSEMBLY ~;~
776-279	AR-279	CATALYST HEATER TUBE 10.75"
776-283	AR-283	LANCE, QUARTZ, CHN600,800
776-284	AR-284	""O"" RING, COMB.TUBE FP228
776-286	AR-286	COMB TUBE CHN600
776-685	AR-6685	TYNA-MITE SOLENOID VALVE
776-923	AR-923	PEDESTAL ONLY
776-931	AR-606	PEDESTAL ~;~
776-962	AR-962	HEATING ELEMENT
776-978	AR-978	GLASS SCOOP, MAC400
776-992	AR-6922	SHIELD & COLLAR
777-060	AR-7060	QUARTZ COMB/COLLECT TUBE
777-064	AR-928	PYREX DELIVERY TUBE SC32 ~;~
777-066	AR-7066	QUARTZ LINER SC32/132
777-073	AR-7073	THERMOCOUPLE FOR USE ON SC32/1
778-114	AR-8114	O RING, CAT HEATER, 111
778-116	AR-8116	O RING, BALLAST
778-120	AR-781	GASKET, BALLAST
778-379	AR-8380	COMBUSTION TUBE INSULATOR
778-492	AR-456	THIS IS NO LONGER AVAILABLE!!!
778-596	AR-596	SST FILTER
778-674	AR-674	GRAR-APHITE CRUCIBLE, 2 GRAM CAR-APA
778-674	AR-674P	GRAR-APHITE CRUCIBLE, 2 G ~;~
778-891	AR-8891	QUARTZ CRUCIBLE, LARGE
778-892	AR-8892	QUARTZ COVER, LARGE
780-044	AR-609	BOAT PULLER
780-047	AR-1202	REACTION TUBE FOR DH102
780-093	AR-093	QUARTZ CRUCIBLE TGA
780-094	AR-094	QUARTZ CRUCIBLE COVER TGA
780-379	AR-379	QUARTZ WINDOW
780-831	AR-831	O-RING
780-853	AR-853	THERMISTOR, LEVEL
780-890	AR-7808	LARGE GRAR-APHITE CRUCIBLE

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
780-892	AR-0892	GRAR-APHITE COVERS, 1000/PKG
780-899	AR-899	HONEYCOMB SUPPORT
780-911	AR-910	CATALYST HTR. TUBE 14.5"
780-927	AR-927	PRINTER PAR-APER ROLL, 4 1/2" X 19
780-928	AR-926	PRINTER RIBBON, RA/SERIES/WR11
780-953	AR-922	PEDESTAL CS-125 ~;~
781-098	AR-1098	THERMISTOR HIGH WATER
781-101	AR-101	THERMISTOR LOW WATER
781-165	AR-469	RELAY
781-273	AR-028	ALPHA CRUCIBLE PLATE, 1" H X 7
781-335	AR-1335	QUARTZ BOAT FOR RC UNITS
781-371	AR-371	FAN FOR CS-244
781-424	AR-1424	O-RING
781-424	AR-1442	O-RING
781-585	AR-1585	OVER TEMP THERMOCUPLE, RC612
781-671	AR-671	INK JET PRINTER PAR-APER
781-883	AR-1883	INK JET PRINTER PAR-APER 500
781-902	AR-902	PRINT CARTRIDGE
782-004	AR-928	PYREX DELIVERY TUBE SC32 ~;~
782-010	AR-7060	QUARTZ COMB/COLLECT TUBE
782-010	AR-7069	REPLACED BY AR7060 ~;~
782-038	AR-2038	COMBUSTION TUBE, RC-412
782-059	AR-2059	NICKEL BOAT LINERS, PK OF 10 L
782-079-110	AR-9023	BOAT STOP, COARSE
782-079-115	AR-9110	CERAMIC BOAT STOP HI FLOW
782-079-115	AR-9115	DISCONTINUED
782-079-120	AR-9022	BOAT STOP HALF SC32/132
782-079-120	AR-9120	DISCONTINUED
782-338	AR-2338	INDICATOR TUBE, OXYGEN
782-338-120	AR-786	OMI PURIFIER TUBE FOR OHN
782-428	AR-2428	CATALYST TUBE
782-534	AR-3930	ELECTRODE INSERT REPAIR, CUSTO
782-601	AR-2601	VACUUM CLEANER 240 V
782-601	AR-2602	VACUUM CLEANER 120 V
782-703	AR-2703	OUTER GRAR-APH CRUC, USE WITH AR2
782-703-100	AR-2703	OUTER GRAR-APH CRUC, USE WITH AR2
782-703-100	AR-2703P	OUTER CRUCIBLE, 100/PKG
782-719	AR-782HD	GRAR-APHITE CRUCIBLE, HIGH TEMP,
782-719	AR-782P	REPLACED BY AR782HDP AND AR782
782-720	AR-782	REPLACED BY AR782HD & AR7827
782-720	AR-7827	GRAR-APHITE CRUCIBLE, HIGH TEMP,
782-720	AR-782HD	GRAR-APHITE CRUCIBLE, HIGH TEMP,
782-720-S	AR-7827	GRAR-APHITE CRUCIBLE, HIGH TEMP,
782-720S	AR-720S	HITEMP SPEC GRAR-APHITE CRUCIBLE
782-720SHD	AR-782HD	GRAR-APHITE CRUCIBLE, HIGH TEMP,
782-720SHD	AR-782P	REPLACED BY AR782HDP AND AR782
782-720SHD	AR-872HD	GrAR-APhite Crucible, High Densit

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
782-721	AR-721	HI-TEMP ELECTRODE INSERT USED
782-721	AR-721LTV	MODIFIED HIGH-TEMP ELECT- ~;~
782-795	AR-2795	INNER GRAR-APH CRUC, USE WITH AR2
782-835	AR-3930	ELECTRODE INSERT REPAIR, CUSTO
782-877	AR-877	TRIM STEP BRUSH 782-918
782-880	AR-2880	VACUUM FILTER, FOR CS UNITS
782-918	AR-877H	HANDLE ONLY FOR AR 877
782-938	AR-938	COMPOSITE BRUSH
782-942	AR-938H	HANDLE ONLY FOR AR 938
783-177	AR-962	HEATING ELEMENT
783-189	AR-3189	"O" RING, CHN1000 (BROWN VI
783-190	AR-3897	O RING, CHN1000, S144
783-193	AR-8343	HEATING ELEMENT, CHN-1000
783-274	AR-278	FURNACE LANCE, CHN1000
783-278	AR-278	FURNACE LANCE, CHN1000
783-299	AR-2038	COMBUSTION TUBE, RC-412
783-301	AR-3301	BRUSH END .38DAI
783-302	AR-3302	BRUSH END MODIFIED
783-515	AR-515	O Ring
783-568	AR-AEE2482	LOWER ELECTRODE TIP
783-785-110	AR-786	OMI PURIFIER TUBE FOR OHN
783-785-HAZ	AR-786	OMI PURIFIER TUBE FOR OHN
783-786	AR-786	OMI PURIFIER TUBE FOR OHN
783-897	AR-3897	O RING, CHN1000, S144
789-094	AR-618	DISCONTINUED
789-095	AR-9095	Spacer for Reagent Tube
789-114	AR-6317	PURGE TUBE, WITH FITTINGS
789-115	AR-789	LANCE TUBE ASSEMBLY, 4.28 INCH
789-116	AR-4174	PURGE TUBE, FOR USE ON SC DETE
789-118	AR-9118	LANCE TUBE ASSEMBLY NO AUTOLOA
789-122	AR-2601	VACUUM CLEANER 240 V
789-412	AR-9412	CATALYST TUBE, 10.75 INCH, TC-
800-635	AR-AM0635	COTTON GLOVES 1 Pair
801-136	AR-AM1136	WHEEL 5" X 0.014" X 1/2"
801-137	AR-AM1137	WHEEL 4" X 0.012" X 1/2", DIAM
801-724	AR-AM1724	Dressing Stick For Diamond Saw
805-481	AR-481	O-RING
808-751	AR-751	O-RING
810-001	AR-AM0001	SiC BELT 4X36 50GRIT W/D 10 Pa
810-002	AR-AM0002	SIC BELT 4X36 80GRIT W/D
810-003	AR-AM0003	SIC BELT 4X36 120GRIT W/D
810-004	AR-AM0004	SIC BELT 4X36 180GRIT W/D
810-005	AR-AM0005	SIC BELT 4X36 240GRIT W/D
810-006	AR-AM0006	SIC BELT 4X36 320GRIT W/D
810-007	AR-AM0007	SIC BELT 4 X 36 400 GRIT W/D
810-008	AR-AM0008	SIC BELT 4 X 36 600 GRIT W/D
810-021	AR-AM0021	ALOX BELT 4 X 36 50 GRIT W/D

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
810-021	AR-AM0022	ALOX BELT 4 X 36 80 GRIT W/D
810-023	AR-AM0023	ALOX BELT 4 X 36 120 GRIT W/D
810-024	AR-AM0024	BELT 4 X 36 180 GRIT W/D
810-025	AR-AM0025	BELT 4 X 36, 240 GRIT W/D
810-033	AR-AM0033	SIC DISC 8 XNH PSA 800 GRIT
810-034	AR-AM0034	SIC DISC 8 XNH PSA 1200 GRIT
810-035	AR-AM0035	8"x NH SiC 800 GRIT 100 Count,
810-036	AR-AM0036	8"x NH SiC 1200 GRIT 100 Count
810-085	AR-AM0085	SiC ABRASIVE ROLL 3"x75' 60 GR
810-086	AR-AM0086	SiC ABRASIVE ROLL 3"x75' 120 G
810-087	AR-AM0087	SiC ABRASIVE ROLL 3"x75' 180 G
810-088	AR-AM0088	SiC ABRASIVE ROLL 3"x75' 240 G
810-089	AR-AM0089	SiC ABRASIVE ROLL 3"x75' 320 G
810-090	AR-AM0090	SiC ABRASIVE ROLL 3"x75' 400 G
810-091	AR-AM0091	SiC ABRASIVE ROLL 3"x75' 600 G
810-103	AR-AM0403	8" Diamond Grinding Disc Plain
810-137	AR-AM1137	WHEEL 4" X 0.012" X 1/2", DIAM
810-161	AR-AM0161	SiC 240 GRIT STRIP 100 Pack
810-162	AR-AM0162	SiC 320 GRIT STRIP 100 Pack
810-163	AR-AM0163	SiC 400 GRIT STRIP 100 Pack
810-164	AR-AM0164	SiC 600 GRIT STRIP 100 Pack
810-219	AR-AM0219	8"x NH SiC 120 GRIT 100 Count
810-221	AR-AM0221	8"x NH SiC 180 GRIT 100 Count
810-223	AR-AM0223	8"x NH SiC 240 GRIT 100 Count
810-223	AR-AM0233	8X.5 SiC, 120 GRIT, 100 Count,
810-225	AR-AM0225	8"x NH SiC 320 GRIT, 100/PK, P
810-227	AR-AM0227	SiC 8XNH PB 400 GRIT
810-229	AR-AM0229	8"x NH SiC 600 GRIT 100 Count,
810-231	AR-AM0231	8X.5", SiC, 60 GRIT, 100 Count
810-235	AR-AM0235	SIC DISC 8 X .5 PB 180 GRIT
810-237	AR-AM0237	SIC DISC 8 X .5 PB 240 GRIT
810-239	AR-AM0239	SIC DISC 8 X .5 PB 320 GRIT
810-241	AR-AM0241	SIC DISC 8 X .5 PB 400 GRIT
810-243	AR-AM0243	SIC DISC 8 X .5 PB 600 GRIT
810-245	AR-AM0245	SIC DISC 12 XNH PSA, 60 GRIT
810-247	AR-AM0247	SIC DISC 12 XNH PSA 120 GRIT
810-249	AR-AM0249	SIC DISC 12 XNH PSA, 180 GRIT,
810-251	AR-AM0251	SIC DISC 12 XNH PSA, 240 GRIT
810-253	AR-AM0253	SIC DISC 12 XNH PSA, 320 GRIT
810-255	AR-AM0255	SIC DISC 12 XNH PSA, 400 GRIT
810-257	AR-AM0256	SIC DISC 12 XNH PSA, 600 GRIT,
810-259	AR-AM0259	SIC DISC 8 XNH PSA, 120 GRIT
810-261	AR-AM0261	SIC DISC 8 XNH PSA, 180 GRIT
810-263	AR-AM0263	SIC DISC 8 XNH PSA, 240 GRIT
810-265	AR-AM0265	SIC DISC 8 XNH PSA, 320 GRIT
810-267	AR-AM0267	SIC DISC 8XNH PSA 400 GRIT
810-269	AR-AM0269	SIC DISC 8 XNH PSA 600GRIT

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
810-271	AR-AM0271	SIC DISC 8 X 2.5 PSA, 120 GRIT
810-273	AR-AM0273	SIC DISC 8 X 2.5 PSA, 180 GRIT
810-275	AR-AM0275	SIC DISC 8 X 2.5 PSA, 240 GRIT
810-277	AR-AM0277	SIC DISC 8 X 2.5 PSA, 320 GRIT
810-279	AR-AM0279	SIC DISC 8 X 2.5 PSA, 400 GRIT
810-281	AR-AM0281	SIC DISC 8 X 2.5 PSA, 600 GRIT
810-291	AR-AM2100	SIC 12XNH PB, 60 GRIT ~;-
810-292	AR-AM2100	SIC 12XNH PB, 60 GRIT ~;-
810-293	AR-AM2103	SIC DISC, 12 X NH, 180 GRIT, P
810-294	AR-AM2104	SIC 12XNH PB 240 GRIT
810-295	AR-AM2105	SIC 12 XNH, PB, 320 GRIT
810-296	AR-AM2106	SIC 12XNH PB 400 GRIT
810-297	AR-AM2107	SIC 12XNH PB 600 GRIT
810-315	AR-AM0467	14" Red Felt Plain Back 10 Pac
810-316	AR-AM0316	14" PAN W PSA 10 Pack
810-327	AR-AM0327	14" Alpha Cloth PSA 10 Pack
810-328	AR-AM0328	14" Red Felt PSA 10 Pack
810-330	AR-AM0333	NYLON 14" " PSA 10 Pack
810-332	AR-AM0332	NYLON 14" " PLAIN BACK 10 Pack
810-333	AR-AM0517	14" Silk PB 10 Pack
810-334	AR-AM0334	14" Alpha Cloth Plain Back 10
810-396	AR-AM0396	8" Diamond Grinding Disc PSA
810-397	AR-AM0397	8" Diamond Grinding Disc PSA
810-398	AR-AM0398	8" Diamond Grinding Disc PSA
810-399	AR-AM0399	8" Diamond Grinding Disc PSA
810-401	AR-AM0401	8" Diamond Grinding Disc PSA
810-402	AR-AM0402	8" Diamond Grinding Disc Plain
810-404	AR-AM0404	8" Diamond Grinding Disc Plain
810-405	AR-AM0405	8" Diamond Grinding Disc Plain
810-406	AR-AM0406	8" Diamond Grinding Disc Plain
810-425	AR-AM0426	14" x NH PSA SiC 180 Grit
810-426	AR-AM0427	14" x NH PSA SiC 240 Grit
810-427	AR-AM0428	14" x NH PSA SiC 320 Grit
810-428	AR-AM0429	14" x NH PSA SiC 400 Grit
810-429	AR-AM0430	14" x NH PSA SiC 600 Grit
810-430	AR-AM0436	ALOX BELT 4X36 50GRIT DRY ONLY
810-437	AR-AM0437	ALOX BELT 4X36 80GRIT DRY ONLY
810-446	AR-AM0425	14" x NH PSA SiC 120 Grit
810-461	AR-AM0461	8" Selvyt Plain Back 10 Pack
810-461	AR-AM0468	8" Velvet Plain Back 10 Pack
810-462	AR-AM0462	NYLON FOR 8" PB 10 Pack
810-463	AR-AM0471	8" Cotton PLAIN BACK 10 Pack
810-465	AR-AM0465	8" Billiard Cloth Plain Back10
810-466	AR-AM0466	8" Red Felt Plain Back 10 Pack
810-467	AR-AM0472	8" Silk PB 10 Pack
810-469	AR-AM0469	CANVAS FOR 8" WHEEL PB 10 Pack
810-470	AR-AM0549	8" Alpha Cloth Plain Back 10 P

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
810-482	AR-AM0482	8" PEDRO FILM W PSA 10 Pack
810-485	AR-AM0485	METAL SAMPLE CLIPS 100 Count
810-489	AR-AM0547	ZIR AL BELT 4X36 180 GRIT
810-491	AR-AM0496	6.5" Alpha Cloth PSA 10 Pack
810-493	AR-AM0548	ZIR AL BELT 4X36 220 GRIT
810-495	AR-AM0495	8.5" Pan W PSA 10 Pack
810-497	AR-AM0497	8" Alpha LAR-AP PSA 10 Pack
810-498	AR-AM0538	8 x 2.5" Alpha Cloth PSA 10 Pa
810-499	AR-AM0546	ZIR AL BELT 4X36 120GRIT
810-510	AR-AM0510	12" Selvyt Plain Back 10 Pack
810-511	AR-AM0511	NYLON FOR 12" WHEEL PB 10 Pack
810-512	AR-AM0512	12" Cotton PLAIN BACK 10 Pack
810-514	AR-AM0514	12" Billiard Cloth Plain Back1
810-515	AR-AM0460	12" Red Felt Plain Back 10 Pac
810-516	AR-AM0516	12" Silk PB 10 Pack
810-518	AR-AM0518	12" CANVAS Plain Back 10 Pack
810-520	AR-AM0539	12" Alpha Cloth Plain Back 10
810-529	AR-AM0529	12" PEDRO FILM W PSA 10 Pack
810-543	AR-AM0540	12" Alpha LAR-AP PSA 10 Pack
810-545	AR-AM0445	8" Imperial Cloth PSA 10 Pack
810-545-010	AR-AM0445	8" Imperial Cloth PSA 10 Pack
810-546	AR-AM0446	12" Imperial Cloth PSA 10 Pack
810-547	AR-AM0451	NYLON 12" " PSA 10 Pack
810-644	AR-AM0644	12" Red Felt PSA 10 Pack
810-645	AR-AM0451	NYLON 12" " PSA 10 Pack
810-647	AR-AM0647	8" Red Felt PSA 10 Pack
810-648	AR-AM0463	NYLON FOR 8" WHEEL PSA 10 Pack
810-665	AR-AM0665	0.05µ Gamma Alumina Suspension
810-697	AR-AM1110	Ultrasonic Cleaner Solution (1
810-778	AR-AM0778	0.05µ Gamma B Aluminum Oxide
810-779	AR-AM0779	0.05µ Gamma B Aluminum Oxide
810-782	AR-AM0782	0.3µ Alpha A Aluminum Oxide 1
810-783	AR-AM0783	0.3µ Alpha A Aluminum Oxide 5
810-786	AR-AM0786	1µ Alpha C Aluminum Oxide 1 Lb
810-787	AR-AM0787	1µ Alpha C Aluminum Oxide 5 Lb
810-788	AR-AM0788	NO CURRENT SOURCE FOR THIS PRO
810-815	AR-AM0486	PLASTIC SAMPLE CLIPS 100 Count
810-818	AR-AM0818	1µ Alumina Suspension 1 Gallon
810-819	AR-AM0832	1µ Alumina Suspension 6 OZ.
810-825	AR-AM0825	5µ Alumina Suspension 6 OZ.
810-825	AR-AM0825G	5µ Alumina Suspension 1 Gallon
810-826	AR-AM0826	0.05µ Gamma Alumina Suspension
810-827	AR-AM0827	0.3µ Alumina Suspension 6 OZ.
810-828	AR-AM0828	DISCONTINUED 3.5.13
810-829	AR-AM0829	DISCONTINUED 3/5/13
810-830	AR-AM0832	1µ Alumina Suspension 6 OZ.
810-831	AR-AM0831	1.0µ Cerium Oxide Suspension 6

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
810-832	AR-AM0834	0.3µ Alumina Suspension 1 Gall
810-847	AR-AM0545	ZIR AL BELT 4X36 60 GRIT W/D
810-856	AR-AM2108	12XNH, SIC, 800 GRIT, PLAIN BA
810-857	AR-AM2109	SIC DISCS, 12XNH, PLAIN PACK,
810-858	AR-AM0858	SIC DISC 12XNH PSA 800 GRIT
810-859	AR-AM0859	SIC DISC 12XNH PSA 1200 GRIT
810-865	AR-AM0865	0.25µ Diamond Compound 5 Grams
810-866	AR-AM0866	0.25µ Diamond Compound 20 Gram
810-867	AR-AM0867	0.625µ Diamond Compound 5 Gram
810-868	AR-AM0868	0.625µ Diamond Compound 20 Gra
810-869	AR-AM0887	1.25µ Diamond Compound 5 Grams
810-870	AR-AM0888	1.25µ Diamond Compound 20 Gram
810-871	AR-AM0889	3µ Diamond Compound 5 Grams
810-872	AR-AM0890	3µ Diamond Compound 20 Grams
810-873	AR-AM0891	6µ Diamond Compound 5 Grams
810-874	AR-AM0892	6µ Diamond Compound 20 Grams
810-880	AR-AM0898	30µ Diamond Compound 20 Grams
810-881	AR-AM0899	45µ Diamond Compound 5 Grams
810-882	AR-AM0901	45µ Diamond Compound 20 Grams
810-885	AR-AM0885	14" x NH PSA SiC 800 Grit
810-886	AR-AM0886	14" x NH PSA SiC 1200 Grit
810-912	AR-AM0893	10µ Diamond Compound 5 Grams
810-913	AR-AM0894	10µ Diamond Compound 20 Grams
810-914	AR-AM0895	15µ Diamond Compound 5 Grams
810-915	AR-AM0896	15µ Diamond Compound 20 Grams
810-938	AR-AM0508	45µ Diamond Spray 4 OZ.
810-939	AR-AM0507	30µ Diamond Spray 4 OZ.
810-940	AR-AM0506	15µ Diamond Spray 4 OZ.
810-941	AR-AM0505	9µ Diamond Spray 4 OZ.
810-942	AR-AM0504	6µ Diamond Spray 4 OZ.
810-943	AR-AM0503	3µ Diamond Spray 4 OZ.
810-944	AR-AM0502	1µ Diamond Spray 4 OZ.
810-945	AR-AM0500	0.25µ Diamond Spray 4 OZ.
810-990	AR-AM3018	1" SILICONE RUBBER CUPS Each
810-991	AR-AM3019	1.25" SILICONE RUBBER CUPS Eac
810-992	AR-AM3020	1.5" SILICONE RUBBER CUPS Each
810-994	AR-AM1018	15µ Diamond Slurry 400 Grams
810-995	AR-AM1016	10µ Diamond Slurry 400 Grams
810-996	AR-AM1014	6µ Diamond Slurry 400 Grams
810-997	AR-AM1012	3µ Diamond Slurry 400 ML
810-998	AR-AM1010	1µ Diamond Slurry 400 Grams
811-002	AR-AM1002	Diamond Lubricant WS 1 Pint
811-002	AR-AM1006	Diamond Lubricant OS 1 Pint
811-003	AR-AM1003	Diamond Lubricant WS 1 Quart
811-003	AR-AM1007	Diamond Lubricant OS 1 Quart
811-004	AR-AM1004	Diamond Lubricant WS 1 Gallon
811-004	AR-AM1008	Diamond Lubricant OS 1 Gallon

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
811-023	AR-AM1023	Cutting Fluid For Diamond Saws
811-024	AR-AM1024	CUTTING FLUID FOR DIAMOND SAWS
811-035	AR-AM1035	WHEEL 9 X 1.25 AL203 10 Pack
811-037	AR-AM1036	WHEEL 9 X1.25 AL203
811-037	AR-AM1142	Blue Diallyl Phthalate (5 Poun
811-043	AR-AM1037	WHEEL 9 X 1.25 AL203 10 Pack
811-049	AR-AM1049	WHEEL 9 X 1.25 AL203
811-053	AR-AM1052	Wheel 9" x 1.25 SiC 10 Pack
811-053	AR-AM1053	WHEEL 9 X 1.25 SiC 10 Pack
811-059	AR-AM1059	WHEEL 10 X 1.25 AL203 10 Pack
811-059	AR-AM1061	WHEEL 10 X 1.25 AL203 10 Pack
811-061	AR-AM1050	WHEEL 10 X 1.25 SIC 10 Pack
811-065	AR-AM1065	WHEEL 12 X 1.25 AL203 10 Pack
811-067	AR-AM1067	WHEEL 10 X1.25 AL203 10 Pack
811-073	AR-AM1051	WHEEL 10 X 1.25 AL203 10 Pack
811-081	AR-AM1081	WHEEL 12 X 1.25 AL203 10 Pack
811-083	AR-AM1083	WHEEL 12 X 1.25 AL2O3 10 Pack
811-089	AR-AM1089	WHEEL 12 X 1.25 AL2O3 10 Pack
811-095	AR-AM1095	WHEEL 12 X 1.25 AL2O3 10 Pack
811-102	AR-AM1066	WHEEL 12 X 1.25 SiC 10 Pack
811-104	AR-AM1104	WHEEL 12 X 1.25 SiC 10 Pack
811-108	AR-AM1108	NITRITE FREE COOLANT (1 Quart)
811-109	AR-AM1109	NITRITE FREE COOLANT(1 Gallon)
811-1111	AR-AM1111	Black Phenolic (5 Pounds)
811-118	AR-AM1118	Red Phenolic (5 Pounds)
811-123	AR-AM1123	Green Phenolic (5 Pounds)
811-128	AR-AM1128	Black Phenolic, Mottled (5 Pou
811-133	AR-AM1133	Transparent Lucite (5 Pounds)
811-135	AR-AM1135	Blue Diallyl Phthalate (5 Poun
811-146	AR-AM1146	NO LONGER AVAILABLE
811-201	AR-AM1201	KOLDMOUNT KIT, 1 POUND RESIN/1
811-202	AR-AM1202	KOLDMOUNT POWDER RESIN 1 POUND
811-203	AR-AM1203	KOLDMOUNT LIQUID HARDENER 12OZ
811-271	AR-AM1271	SILICONE MOLD RELEASE 16 Oz. H
811-275	AR-AM1275	GRADUATED PAR-APER CUPS, 5 OZ
811-279	AR-AM1279	GRADUATED PAR-APER CUPS 10 OZ.
811-301	AR-AM1301	BLACK 1" PREFORMS 100 Count
811-302	AR-AM1302	BLACK 1" PREFORMS 500 Count
811-303	AR-AM1303	BLACK 1" PREFORMS 1000 Count
811-305	AR-AM1305	BLK 1 1/4" Preforms 100 Count
811-306	AR-AM1306	BLK 1 1/4" PREFORMS 500 Count
811-307	AR-AM1307	BLK 1 1/4" PREFORMS 1000 Count
811-309	AR-AM1309	BLK 1 1/2" PREFORMS 100 Count
811-310	AR-AM1310	BLK 1 1/2" PREFORMS 500 Count
811-311	AR-AM1311	BLK 1 1/2" PREFORMS 1000 Count
811-313	AR-AM1313	RED 1" PREFORMS 100 Count
811-314	AR-AM1314	RED 1" PREFORMS 500 Count

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
811-315	AR-AM1315	RED 1" PREFORMS 1000 Count
811-317	AR-AM1317	RED 1 1/4" PREFORMS 100 Count
811-318	AR-AM1318	RED 1 1/4" PREFORMS 500 Count
811-319	AR-AM1319	RED 1 1/4" PREFORMS 1000 Count
811-321	AR-AM1321	RED 1 1/2" PREFORMS 100 Count
811-322	AR-AM1322	RED 1 1/2" PREFORMS 500 Count
811-323	AR-AM1323	RED 1 1/2" PREFORMS 1000 Count
811-325	AR-AM1325	GREEN 1" PREFORMS 100 Count
811-326	AR-AM1326	GREEN 1" PREFORMS 500 Count
811-327	AR-AM1327	GREEN 1" PREFORMS 1000 Count
811-329	AR-AM1329	GREEN 1 1/4" PREFORMS 100 Coun
811-330	AR-AM1330	GREEN 1 1/4" PREFORMS 500 Coun
811-331	AR-AM1331	GREEN 1 1/4" PREFORMS 1000 Cou
811-333	AR-AM1333	GREEN 1 1/2" PREFORMS 100 Coun
811-334	AR-AM1334	GREEN 1 1/2" PREFORMS 500 Coun
811-335	AR-AM1335	GREEN 1 1/2" PREFORMS 1000 Cou
811-563	AR-AM1161	EPOXY KIT HMA HMG
811-563	AR-AM1162	EPOXY KIT HMA HMG
811-563	AR-AM1163	EPOXY RESIN, 1 Quart
811-563	AR-AM1164	EPOXY RESIN, 1 Gallon
811-563	AR-AM1165	EPOXY HARDENER, 8 OZ, HMA, HM
811-563	AR-AM1167	FAST SET EPOXY RESIN, 1 QT
811-563	AR-AM1168	FAST SET EPOXY HARDENER HMA H
811-563	AR-AM1169	FAST SET EPOXY RESIN, 1 GALLON
811-563	AR-AM1170	FAST SET EPOXY HARDENER, 1 QT,
811-634	AR-AM1634	8" AlphaGrab Disc 10 Pack
811-635	AR-AM1635	12" AlphaGrab Disc 10 Pack
812-021	AR-AM0341	8" x NH PSA Zir AI 60 Grit
812-024	AR-AM0342	8" x NH PSA Zir AI 120 Grit
812-045	AR-AM0357	12" x NH PSA Zir AI 60 Grit
812-046	AR-AM0358	12" x NH PSA Zir AI 120 Grit
812-110	AR-AM1022	0.06µ Master Aqueous Polish 32
812-121	AR-AM1022S	0.06µ Master Aqueous Polish 1
812-201	AR-AM2120	SIC DISCS, 10XNH, PLAIN BACK,
812-202	AR-AM2090	SIC DISC, 10 X NH, PSA, 60 GRIT
812-203	AR-AM2122	SILICON CARBIDE, 120 GRIT
812-204	AR-AM2091	SIC DISCS, 10XNH, PSA, 120 GRI
812-205	AR-AM2123	SIC DISCS, 10XNH, PLAIN BACK,
812-206	AR-AM2092	SIC DISC 10 X NH PSA 180 GRIT
812-207	AR-AM2124	SIC DISC, 10XNH, PLAIN BACK, 2
812-208	AR-AM2093	SIC DISCS, 10XNH, 240GRIT, PSA
812-209	AR-AM2125	SIC DISCS, 10XNH, PLAIN BACK,
812-211	AR-AM2126	SiC DISC, 10XNH PB, 400 GRIT
812-212	AR-AM2095	SIC DISCS, 10XNH, 400GRIT, PSA
812-213	AR-AM2127	SIC DISC 10 X NH PB 600 GRIT
812-214	AR-AM2096	SIC DISC, 10XNH, 600 GRIT, PSA
812-217	AR-AM2128	SIC DISCS, 10XNH, PLAIN BACK,

Cross references: OEM & Alpha Resources Article Numbers

LECO® Article No.	Alpha Resources Article No.	Description
812-218	AR-AM2097	SIC DISCS, 10XNH, PSA BACK, 80
812-219	AR-AM2129	SIC DISCS, 10XNH, PLAIN BACK,
812-220	AR-AM2098	SIC DISCS, 10XNH, PSA BACK, 12
813-090	AR-AM0995	STIRRING STICKS 500 Count
820-200	AR-AM0200	Ultrasonic Cleaner
874-000-181	AR-4181	Small Nickel Boat
906-116-110	AR-9116	COMB TUBE NEW STYLE, FP 428, W
UK-1992	AR-AEC4050	OXYGEN LANCE STRAIGHT TRANSPAR
UK-7750	AR-AEC4048	CRUCIBLE TRANSPARENT SILICA VE