

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
4045	AR-04045	ELTRA TRANSFORMER
5000	AR-05000	ELTRA IR CELL
5020	AR-05020	ELTRA CHOPPER MOTOR BLADE
5030	AR-05030	ELTRA CHOPPER BLADE HOLDER
5048	AR-05048	ELTRA INFRARED SOURCE (EMITTER
5060	AR-05060	ELTRA REFLECTOR
5064	AR-05064	ELTRA CHOPPER MOTOR
5065	AR-05065	ELTRA CHOPPER MOTOR
5150	AR-05150	ELTRA IR PATH CONNECTOR WITH W
5160	AR-05160	ELTRA IR PATH CONNECTOR WITH W
5170	AR-05170	ELTRA GAS CONDITIONER
5244	AR-05244	ELTRA IR 4 CHOPPER
5251	AR-05251	CHOPPER MODULE WITH 3 IR SOURC
5260	AR-05260	DO NOT SELL - ELTRA IR PATH TU
5270	AR-05270	ELTRA THREADED ROD
5275	AR-05275	**OBSOLETE** USE AR05274
05502-295	AR-05502-295	ELTRA IR PATH
05527-1	AR-05527-1	ELTRA IR PATH
05542-10	AR-05542-100	ELTRA IR PATH
05547-20	AR-05547-20	ELTRA IR PATH
6005	AR-06005	ELTRA IR 21 INFRARED ELECTRONI
6025	AR-06025	ELTRA IRT ZERO ADJUSTMENT BOARD
6032	AR-06032	ELTRA INFRARED DETECTOR WITH C
6034	AR-06034	IR DETECTOR W/ SO2 FILTER
6040	AR-06040	REPLACED BY AR06421, IR POWER
6058	AR-06058	ELTRA TEMPERATURE CONTROL BOAR
6059	AR-06059	ELTRA REAR PANEL WITH COMPONENT
6210	AR-06210	ELTRA IR CABLE
6251	AR-06251	ELTRA IR CELL TRANSFORMER
6421	AR-06421	ELTRA INFRARED POWER SUPPLY AS
6670	AR-06670	ELTRA IR ZERO ADJUSTMENT BOARD
6714	AR-06741	ELTRA INFRARED ELECTRONICS BOA
9090	AR-09090	ELTRA REAGENT TUBE
9271	AR-09271	REPLACED BY AR15270
10169	AR-10169	ELTRA GLASS SUPPORT
10196	AR-10196	ELTRA GLASS SUPPORT
11035	AR-11035	ELTRA COOLING FAN
11042	AR-11042	ELTRA UPPER TUBE HOLDER
11045	AR-11045	ELTRA LOWER TUBE HOLDER
11062	AR-11062	ELTRA REAGENT TUBE 160MM X 16M
11064	AR-11064	REPLACED BY AR88400-0006
11064	AR-88400-0006	ELTRA REAGENT TUBE 280MM X 16M

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
11091	AR-11091	ELTRA DUST TRAP MECHANISM
11093	AR-11093	ELTRA DUST FILTER PLUG
11105	AR-311	SST FILTER
11106	AR-11106	ELTRA DUST FILTER CARTRIDGE
11110	AR-11110	ELTRA FILTER HOUSING FOR AR 31
11115	AR-73040	CONNECTION TUBE 6X9MM 1 METER
11120	AR-11120	ELTRA PAPER FILTER HOLDER
11170	AR-980	IN LINE PARTICLE FILTER
11180	AR-980	IN LINE PARTICLE FILTER
11185	AR-11185	ELTRA PAPER FILTERS, PKG OF 10
11190	AR-11190	ELTRA EXHAUST MUFFLER
11370	AR-11370	SEE AR11000-2001 FOR CS800 AND
11376	AR-11376	ELTRA FIXING PLATE FOR AR11822
11380	AR-11380	ELTRA PNEUMATIC VALVE FOR FURN
11390	AR-11390	ELTRA OXYGEN SOLENOID VALVE
11400	AR-11400	ELTRA PRESSURE OUTLET SOLENOID
11408	AR-11408	ELTRA PRESSURE SWITCH
11410	AR-11410	OXYGEN STOP SOLENOID VALVE
11415	AR-11415	ELTRA OXYGEN STOP SOLENOID VAL
11420	AR-11420	ELTRA COIL COOLING SOLENOID VA
11430	AR-11430	ELTRA VALVE V3
11440	AR-11440	ELTRA BYPASS SOLENOID VALVE
11470	AR-11470	ONE-WAY FILTER VALVE
11480	AR-11480	ELTRA ADJUSTABLE RESTRICTOR
11492	AR-11492	ELTRA PRESSURE REGULATOR
11495	AR-11495	ELTRA PURGE PRESSURE REGULATOR
11822	AR-11822	ELTRA PNEUMATIC VALVE BLOCK
12016	AR-12016	ELTRA GAS FLOW AND FURNACE CON
12044	AR-12044	ELTRA TRANSFORMER
12045	AR-12045	ELTRA TRANSFORMER
12080	AR-12080	ELTRA RECTIFIER
12100	AR-12100	ELTRA TRANSFORMER
13060	AR-13060	ELTRA PROTECTION SHEET
13067	AR-13067	ELTRA COMBUSTION COIL
13080	AR-13080	CAPACITOR SUPPORT
13081	AR-13081	ELTRA CAPACITOR SUPPORT
13090	AR-13090	ELTRA UPPER COIL CONNECTOR
13100	AR-13100	ELTRA LOWER COIL CONNECTOR
13110	AR-13110	ELTRA ANODE CONNECTOR
13120	AR-13120	TANK STRAP
13122	AR-13122	ELTRA CAPACITOR CONNECTOR
13130	AR-13130	ELTRA CAPACITOR CONNECTOR

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
13150	AR-13150	ELTRA ANODE HEAT SINK
13160	AR-13160	ELTRA COIL HEAT SINK
13170	AR-13170	ELTRA RADIATION SHIELD, REPLAC
13175	AR-13175	ELTRA INSULATOR
13210	AR-13210	ELTRA GRID CHOKE
13220	AR-13220	ELTRA ANODE CHOKE
13250	AR-13250	ELTRA CHASSIS SUPPORT
13260	AR-13260	ELTRA HIGH VOLTAGE FILTER
13270	AR-13270	ELTRA RESISTOR
13850	AR-13850	ELTRA CENTRIFUGAL BLOWER WITH
14009	AR-14009	ELTRA PNEUMATIC CYLINDER FOR F
14014	AR-14014	ELTRA COMPLETE FURNACE CLEANIN
14021	AR-14021	ELTRA UPPER FURNACE LOCK
14025	AR-14025	LOWER FURNACE BLOCK
14026	AR-14026	INCORRECT PART NUMBER - SEE AR
14030	AR-14030	ELTRA UPPER FURNACE PLATE
14035	AR-14035	ELTRA LOWER FURNACE PLATE
14041	AR-14041	LANCE WASHER
14045	AR-14045	ELTRA CLEANING BRUSH FOR COMBU
14050	AR-14050	BRUSH HOLDER
14051	AR-14051	ELTRA BRUSH HOLDER
14072	AR-14072	ELTRA CERAMIC HEAT SHIELD FOR
14080	AR-14080	ELTRA CLEANING MECHANISM ROD
14090	AR-14090	ELTRA BEARING
14100	AR-14100	ELTRA MOUNTING ROD
14110	AR-14110	ELTRA THREADED ROD
14120	AR-14120	ELTRA FURNACE COVER
14130	AR-14130	ELTRA/BRUKER QUARTZ COMBUSTION
14140	AR-14140	LOADING TRAY
14160	AR-14160	ELTRA UPPER LOG
14161	AR-14161	ELTRA LOWER KNURLED NUT, CS800
14165	AR-14165	ELTRA UPPER KNURLED NUT, CS800
14168	AR-14168	ELTRA PEDESTAL
14170	AR-14170	ELTRA PEDESTAL MOUNT
14180	AR-14180	ELTRA FURNACE CLOSURE
14185	AR-14185	ELTRA TRAY
14200	AR-14200	ELTRA METAL TUBE
14210	AR-14210	ELTRA THREADED ROD
14220	AR-14220	ELTRA CYLINDER SUPPORT
14531	AR-14531	AUTOLOADER FOR 104 POSITIONS
14532	AR-14532	DO NOT SELL - SEE AR14530-1001
14535	AR-14535	ELTRA AUTOLOADER FOR HELIOS 36

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
14647	AR-14647	LASER FOR AUTOSAMPLER
15005	AR-15005	GASFLOW CONTROL BOARD
15037	AR-15037	ELTRA VOLTAGE STABILIZER CARD
15038	AR-15038	ELTRA PUMP CONTROL BOARD NR 3
15063	AR-15063	GASFLOWMETER PHOTO 300L/H
15083	AR-15083	ELTRA GAS FLOW INDICATOR 15 I/
15085	AR-15085	ELTRA GAS FLOW INDICTOR 130H/L
15087	AR-15087	ELTRA FLOW DISPLAY 300I/H
15095	AR-15095	ELTRA FLOW DISPLAY 600I/H
15098	AR-15098	ELTRA FC 21 PUMP CONTROL BOARD
15099	AR-15099	ELTRA GAS FLOW CONTROL BOARD F
15115	AR-15115	ELTRA GAS FLOW SENSOR
15116	AR-15116	FLOW SENSOR, FOR CS-800 REPLAC
15150	AR-15150	FLOW CONTROLLER
15155	AR-15155	PRESSURE SENSOR
15170	AR-15170	ELTRA GAS FLOW SENSOR
15175	AR-15175	ELTRA GAS FLOW SENSOR
15207	AR-15207	ELTRA GAS FLOW SENSOR
15209	AR-15209	ELTRA FLOW SENSOR
15211	AR-15211	FLOW SENSOR, CS 800 REPLACES A
15212	AR-15212	FLOW SENSOR CONTROLLER, CS-500
15220	AR-15220	PRESSURE TRANSDUCER ELTRA CS-8
15222	AR-15222	PRESSURE TRANSDUCER ELTRA CS58
15266	AR-15266	PUMP, M2K3
15267	AR-15267	REPLACED BY AR15270
15268	AR-15268	ELTRA PUMP CONTROL BOARD PC
15270	AR-15270	SEE AR88400-0283 - DO NOT SELL
15275	AR-15275	PUMP SUPPORT
16100	AR-16100	ELTRA POWER SUPPLY BOARD NK 31
16110	AR-16110	ELTRA POWER SUPPLY BOARD NK31
16113	AR-16113	ELTRA POWER SUPPLY
16114	AR-16114	ELTRA 24VDC POWER SUPPLY
17020	AR-17020	CTLST. FURN. CONT. BD.NR5
18020	AR-18020	CONTROL KEYBOARD
18150	AR-18150	NO LONGER NEEDED ON ELTRA
18160	AR-18160	BALANCE INTERFACE CABLE
18171	AR-18171	COMPUTER INTERFACE CABLE
18173	AR-18173	ELTRA PC CABLE 9-PIN SERIAL
18175	AR-18175	DISCONTINUED
18175X	AR-18175X	ELECTRONIC RACK EXCHANGE
18454	AR-18454	ELTRA BALANCE CABLE
18467	AR-18467	ELTRA MICROCONTROLLER BOARD 18

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
200000	AR-20000	SEE AR20002
20002	AR-20002	ELTRA CATALYST FURNACE
20036	AR-20036	ELTRA UPPER TUBE HOLDER
20040	AR-20040	ELTRA CATALYST TUBE 200 X 16MM
20045	AR-20045	ELTRA LOWER TUBE HOLDER
21010	AR-21010	ELTRA HEATER SECTION
21022	AR-21022	PURIFYING FURNACE LOCK
21050	AR-21050	PURIFYING FURNACE HEATER
21060	AR-21060	PURIF.FURN. MNTG.PLATE
21070	AR-21070	SMALL PIPE
21080	AR-21080	SMALL ROD
21090	AR-21090	PURIFYING FURNACE PIPE
21100	AR-21100	PURIFYING FURNACE ROD
21110	AR-21110	PURIFYING FURNACE HEATER
21120	AR-21120	ELTRA QUARTZ TUBE
21130	AR-21130	PUR. FURN. SCREEN GUARD
24108	AR-24108	NO LONGER VALID
24250	AR-24250	ELTRA VALVE V10
24260	AR-24260	ELTRA VALVE V13
24265	AR-24265	ELTRA VALVE V12
24285	AR-24285	ELTRA BOARD CS20
24310	AR-24310	VITON TUBE 6X4X1MM
25105	AR-25105	NO LONGER VALID
25186	AR-25186	REPLACED BY AR25188, INLET VAL
25188	AR-25188	INLET VALVE, REPLACES AR25186
25190	AR-25190	ELTRA VALVE V8 230 VAC COIL
25195	AR-25195	SOLENOID VALVE, ONH2000
25630	AR-25630	EXTERNAL GAS PURIFIER FOR ELTR
26063	AR-26063	ELTRA CRUCIBLES FOR THEMOSTEP
26319	AR-26319	ELTRA HOLLOW SCREW
26618	AR-26618	ELTRA VALVE BLOCK
31205	AR-31205	ELTRA UPPER FURNACE PART
31220	AR-31220	ELTRA LOWER FURNACE ASSEMBLY
31230	AR-31230	ELTRA UPPER FURNACE PART
31246	AR-31246	ELTRA COMBUSTION CHAMBER
31250	AR-31250	ELTRA UPPER ELECTRODE INSERT
31322	AR-31322	ELTRA UPPER FURNACE PART
31325	AR-31325	ELTRA SAMPLE DROP MECHANISM
31331	AR-31331	ELTRA COOLING COVER
31345	AR-31345	ELTRA UPPER ELECTRODE INSERT
31350	AR-31350	LOWER ELECTRODE
31355	AR-31355	GRAPHITE INSERT

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
31360	AR-31360	ELTRA GRAPHITE TIP
31365	AR-31365	DO NOT SELL - USE AR31365-8000
31380	AR-31380	ELTRA ISOLATING RING
31393	AR-31393	ELTRA LOWER POWER CONNECTOR
31411	AR-31411	ELTRA UPPER POWER CONNECTOR
31465	AR-31465	ELTRA ECCENTRIC WHEEL
31476	AR-31476	MOTOR ASSEMBLY FOR ON/OH/ONH
31480	AR-31480	ELTRA SUPPORT FOR MOTOR BRACKE
32066	AR-32066	ELTRA PW 13 FURNACE CONTROL BO
32210	AR-32210	TEMP CONTROL BOARD ON 21
32400	AR-32400	ELTRA TRANSFORMER
32450	AR-32450	ELTRA BRAIDED CABLE FOR ONH SE
33075	AR-33075	SEE AR25196
33105	AR-33105	NO LONGER VALID
33400	AR-33400	COOLING TANK
33420	AR-33420	WATER TANK FOR ELTRA ONH UNITS
33440	AR-33440	ELTRA WATER LEVEL TANK
33485	AR-33485	ELTRA WATER TEMPERATURE SENSOR
33500	AR-33500	SEE AR33501-2001
33515	AR-33515	PUMP FITTINGS
33540	AR-33540	DO NOT USE - SEE AR35540
33650	AR-33650	TAP WATER VALVE
35038	AR-35038	CONTROL BOARD ASSEMBLY FOR ELT
35061	AR-35061	ELTRA PUMP SUPPORT
35310	AR-35310	ELTRA PRESSURE SWITCH
35330	AR-35330	ELTRA ATTENUATOR VOLUME
35336	AR-35336	ELTRA ATTENUATOR VOLUME
35355	AR-35355	ELTRA LOWER REAGENT TUBE CONNE
35541	AR-35541	ELTRA CERAMIC LANCE, HELIOS
35542	AR-AEC6512	CERAMIC PEDESTAL
35550	AR-35550	NO LONGER VALID
35782	AR-35782	POWER SUPPLY BOARD
35785	AR-35785	CABLE
36001	AR-36001	ELTRA FURNACE FOR CS-500
36101	AR-36101	ELTRA BOAT STOP
36120	AR-1335	QUARTZ BOAT FOR RC UNITS
36216	AR-36216	ELTRA COMBUSTION BOAT INSERTIO
36217	AR-36217	ELTRA COMBUSTION BOAT REMOVING
36240	AR-36240	ELTRA CERAMIC PLATE
36242	AR-36242	ELTRA CERAMIC PLATE
36283	AR-36283	ELTRA HEATING ELEMENT CONNECTO
36288	AR-36288	REPLACED BY AR36289

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
36289	AR-36289	ELTRA HEATING ELEMENT CONNECTO
36445	AR-36445	ELTRA PLATFORM FOR COMBUSTION
36465	AR-36465	ELTRA CERAN PLATE
36620	AR-36620	ELTRA DUST TRAP
36665	AR-36665	ELTRA SPRING
36700	AR-36700	ELTRA UPPER REAGENT TUBE CONNE
36720	AR-36720	ELTRA OUTLET TUBE FOR RESISTAN
36725	AR-36725	ELTRA TEFLON TUBE
36730	AR-36730	ELTRA TEFLON TUBE 0.8 METER 4m
36735	AR-36735	ELTRA - SEE AR88400-0025 AND O
36751	AR-36751	ELTRA PANEL
36798	AR-36798	ELTRA TRANSFORMER
36799	AR-36799	ELTRA COOLING FAN
36904	AR-36904	ELTRA TEMPERATURE CONTROL BOAR
36907	AR-36907	ELTRA TEMPERATURE CONTROL BOAR
36908	AR-36908	ELTRA CABLE FOR TH 43/TH44
36914	AR-36914	ELTRA THERMOCOUPLE 30CM CABLE
36975	AR-36975	ELTRA FURNACE DISPLAY
36999	AR-36999	ELTRA HTF 540 FURNCE FOR BURNI
38001	AR-38001	ELTRA TIC MODULE FOR USE WITH
38220	AR-38220	GLASS STOPPER, CONNECTION TO T
38225	AR-38225	GLASS STOPPER, CONNECTION TO T
38227	AR-38227	GLASS STOPPER, CONNECTION TO T
38340	AR-38340	GLASS SUPPORT
44500	AR-44500	DISCONTINUED, SEE AR13850
48067	AR-48067	ELTRA LOWER HOLDER FOR REAGENT
48069	AR-48069	ELTRA UPPER HOLDER FOR REAGENT
48070	AR-48070	ELTRA UPPER HOLDER FOR REAGENT
48102	AR-48102	ELTRA METAL DUST FILTER CW800
48105	AR-48105	ELTRA FILTER SUPPORT
48248	AR-48248	ELTRA DUST BOX
48390	AR-48390	ELTRA IR CELL
48750	AR-48750	COMBUSTION TUBE, CW800
48756	AR-48756	RUBBER SLEEVE FOR CW800
48917	AR-48917	ELTRA HEATED TUBE
48960	AR-48960	ELTRA TRANSFORMER
5330	AR-5330	FOUR WAY VALVE
60234	AR-60234	ELTRA PNEUMATIC VALVE
60236	AR-60236	PRESSURE REGULATOR REG-010-1/8
60240	AR-60240	ELTRA SILENCER
60513	AR-60513	ELTRA FITTING 4-1/8"
60522	AR-60522	ELTRA FITTING

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
60525	AR-60525	L FITTING 4-1/4"
60540	AR-60540	I FITTING 4-4
60552	AR-60552	T RESTRICTOR, 4-1/8"
60558	AR-60558	T-FITTING 4-4-1/8"
60566	AR-60566	I FITTING 4-1/4"
60631	AR-60631	ELTRA CLOSURE
60635	AR-60635	ELTRA SEAL FOR AR11822 MANIFOL
66200-0106	AR-66200-0106	ELTRA O-RING
66200-0107	AR-66200-0107	ELTRA O-RING
66200-0108	AR-66200-0108	ELTRA O-RING 51X3 ELEMENTRAC O
66200-0109	AR-66200-0109	O-Ring - Eltra Elementrac
66200-0115	AR-66200-0115	ELTRA O-RING 40X6 ELEMENTRAC O
66200-0159	AR-66200-0159	ELTRA O-RING 12.5X3 ELEMENTRAC
66200-0160	AR-66200-0160	ELTRA O-RING 14.5X3 ELEMENTRAC
66200-120	AR-66200-0120	ELTRA O-RING 56X2 ELEMENTRAC O
70030	AR-70030	SET OF O RINGS FOR ELTRA CS-80
70120	AR-70120	ELTRA O RING FOR CLEANING ASSY
70150	AR-70150	O RING
70180	AR-70180	ELTRA O RING MULTIPIECE IR PAT
70182	AR-2266	FIRING FUSE WIRE, 60 INCH
70210	AR-70210	ELTRA O RING CATALYST FURNACE
70230	AR-70230	ELTRA O RING REAGENT TUBE
70240	AR-70240	ELTRA O-RING
70270	AR-70270	O RING
70280	AR-70280	ELTRA O RING CHOPPER MOTOR
70285	AR-70285	ELTRA O RING
70320	AR-70320	ELTRA O RING INSIDE DUST TRAP
70330	AR-70330	ELTRA O RING IR SOURCE AND PAT
70350	AR-70350	ELTRA O RING DUST TRAP SEAL
70370	AR-70370	ELTRA O RING DUST TRAP FILTER
70380	AR-70380	O RING FOR COMBUSTION TUBE
70390	AR-70390	ELTRA O RING PEDESTAL SEAL
70405	AR-70405	ELTRA O RING ONH LOWER FURNACE
70410	AR-70410	ELTRA O RING ONH SAMPLE DROP A
70415	AR-70415	ELTRA O RING ONH UPPER FURNACE
70425	AR-70425	ELTRA O RING ONH UPPER FURNACE
70427	AR-70427	O RING, HELIOS
70428	AR-70428	O RING, HELIOS
70430	AR-70430	ELTRA SEAL ONH SAMPLE DROP C C
70435	AR-70435	ELTRA SEAL ONH SAMPLE DROP C C
71001	AR-1300	NO LONGER VALID - PLEASE SEE A
71010	AR-71010	ELTRA CLEANING BRUSH FOR PEDES

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
71029	AR-71029	ELTRA FURNACE CLEANING BRUSH F
71031	AR-14072B	ELTRA CLEANING BRUSH
71036	AR-71030	CLEANING BRUSH FOR DUST TRAP
71060	AR-71060	ACID BOTTLE
71061	AR-71061	ACID BOTTLE, 500 ML
71065	AR-71065	DISPENSER ASSIMAT
71070	AR-71070	HEATER WITH MAGNETIC STIRRER
71090	AR-71090	ELTRA VOLTAGE REGULATOR, 5KVA,
72010	AR-72010	ELTRA PRESSURE GAUGE
72020	AR-72020	ELTRA PRESSURE GAUGE
72060	AR-72060	REPLACED BY AR14009, CYLINDER
73020	AR-73020	TRANSPARENT TUBE, 5 METER
73030	AR-73030	ELTRA CONNECTION TUBE 4X6MM 1
73040	AR-73040	CONNECTION TUBE 6X9MM 1 METER
75120	AR-75120	ELTRA SPRING
75122	AR-75122	ELTRA SPRING
75130	AR-75130	ELTRA RETAINING WASHER
75140	AR-75140	ELTRA SAFETY RING
75150	AR-75150	ELTRA METAL TUBE
75190	AR-75190	ELTRA WASHER
76002	AR-76002	ELTRA WASHER, CS800/2000
76005	AR-76005	ELTRA KNURLED NUT, CS800/2000
76010	AR-76010	DUST TRAP KNOB
77030	AR-77030	CURRENT TRANSFORMER
77032	AR-77032	ELTRA CIRCUIT BREAKER 20A
77050	AR-77050	ELTRA TRIAC
77051	AR-77051	ELTRA THYRISTOR BLOCK
77135	AR-77135	ELTRA CAPACITOR
77140	AR-77140	ELTRA HF FILTER, 250V / 16A
77210	AR-354	OSCILLATOR TUBE CS144-600, ELT
77320	AR-77320	CAPACITOR, 100PF/30KVA
77321	AR-77321	CAPACITOR
77330	AR-77330	ELTRA CAPACITOR
77333	AR-77333	ELTRA FILTER
77335	AR-77335	ELTRA CAPACITOR
77337	AR-77337	CAPACITOR
77337	AR-77337	CAPACITOR
77340	AR-77340	REPLACED BY AR77345
77341	AR-77341	CAPACITOR
77350	AR-77350	ELTRA CAPACITOR 100 nF
77410	AR-77410	ELTRA PANEL METER 5A
77423	AR-77423	ELTRA TEMPERATURE SENSOR FOR I

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
77425	AR-77425	ELTRA THERMOCOUPLE 100CM CABLE
77430	AR-77430	ELTRA PANEL METER 100 microA
77450	AR-77450	ELTRA TEMPERATURE SWITCH
77460	AR-77460	ELTRA MICRO SWITCH
77489	AR-77489	ELTRA BRACKET
77490	AR-77490	REPLACED BY AR36289, CLAMP/WIR
77501	AR-77500	ELTRA HEATING ELEMENT SET OF 4
77505	AR-77505	ELTRA SPACER FOR AR77500
77507	AR-77507	REPLACED BY AR77505
77510	AR-77510	ELTRA HEATER FOR INFRA-RED CEL
77525	AR-77525	REPLACED BY 77500, KIT REQUIRE
77600	AR-77600	ELTRA RESISTOR
77610	AR-77610	ELTRA RESISTOR
78010	AR-78010	ELTRA POWER SWITCH - 3 POSITIO
78015	AR-78015	ELTRA MAINS POWER SWITCH 3 POS
78017	AR-78017	ELTRA MAINS POWER SWITCH
78040	AR-78040	ELTRA BUTTON FOR LEAKAGE TEST
78056	AR-78056	MOTOR FOR NEW UNITS - 4 YEARS
80925	AR-80925	SET OF 2 FUSES
80926	AR-80926	SET OF 3 FUSES CS-500
80927	AR-80927	ELTRA SET OF FUSES
82011	AR-82011	FILTER FLEECE (M120)
82120	AR-82120	FAN FILTER AND HOUSING CS800
82130	AR-82130	ELTRA FILTER
88400-0008	AR-008	ALPHACEL II TUNGSTEN/TIN, -20/
88400-0066	AR-AND3500	NICKEL CAPSULE, 12.5 X 5MM, 65
88400-0175	AR-182	EVAC PIN TUBE 150MM X 4MM ID,
88400-0176	AR-8030	FILTERING CRUCIBLE, LARGE, 100
88400-0221	AR-303	NODULAR IRON C=3.86% S=0.0210%
88400-0283	AR-15270	SEE AR88400-0283 - DO NOT SELL
88400-0329	AR-321	ALPHACAT 250g PER BOTTLE
88400-0329	AR-321K	ALPHACAT 1KG PACKAGE
88400-0474	AR-2186-200	TIN FOIL CUPS, 36mm x 36mm
88400-0535	AR-912	PLATINIZED SILICA
88400-0912	AR-027	ALPHACEL TUNGSTEN III -20+40 M
88600-0008	AR-426	COMBUSTION CONTROL, SOLID, I P
88600-0010	AR-173	ALPHACEL II HP, 750G BOTTLE
88600-0010-0005	AR-008	ALPHACEL II TUNGSTEN/TIN, -20/
88600-0011	AR-8053	COMBUSTION BOAT, BOX OF 500
88600-0012	AR-2344	NI BASKETS LOW OXYGEN, 1g, 10
88600-0013	AR-673	WROUGHT IRON CHIP STANDARD C=0
88600-0014	AR-9880S	10 MM HOLE COVER, BAG OF 250

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
88600-0015	AR-1265	CRUCIBLE HOLDER TRAY
88600-0017	AR-040	COVER 4MM HOLE 1000/BOX
90145	AR-929	CRUCIBLE TONGS FOR AR3818
90148	AR-3818	CERAMIC CRUCIBLE, BAGGED, BOX
90148	AR-3818F	CERAMIC CRUCIBLE, BOX OF 1000,
90149	AR-3818F	CERAMIC CRUCIBLE, BOX OF 1000,
90150	AR-3818F	CERAMIC CRUCIBLE, BOX OF 1000,
90152	AR-3818	CERAMIC CRUCIBLE, BAGGED, BOX
90152	AR-3818F	CERAMIC CRUCIBLE, BOX OF 1000,
90153	AR-8203	COMBUSTION BOAT, BOX OF 500
90155	AR-8203	COMBUSTION BOAT, BOX OF 500
90160	AR-8053	COMBUSTION BOAT, BOX OF 500
90162	AR-90162	ELTRA COMBUSTION TUBE
90180	AR-431P	INNER CARBON CRUCIBLE, 100 PER
90180	AR-5431	INNER GRAPHITE CRUCIBLE, BOX O
90180	AR-5431S	INNER GRAPHITE CRUCIBLE THICK
90185	AR-433	OUTER CARBON CRUCIBLE 50 PER P
90185	AR-5433	see AR433 OUTER CARBON CRUCIBL
90190	AR-053HD	GRAPHITE CRUCIBLE, 1000/PACK
90190	AR-053HD-400	GRAPHITE CRUCIBLES 400 PACK
90190	AR-053P	GRAPHITE CRUCIBLE, BAG OF 200
90190	AR-6247	GRAPHITE CRUCIBLE, 1000/PACK
90191	AR-90191	ELTRA GRAPHITE CRUCIBLE
90192	AR-331S	GRAPHITE CRUCIBLE .625 OD, SOL
90200	AR-171	ALPHADRI MAGNESIUM PERCHLORATE
90210	AR-2174	ALPHASOLVE II, 20-30 MESH 500
90210	AR-2176	ALPHASOLVE II, 8-20 MESH 500 G
90220	AR-027	ALPHACEL TUNGSTEN III -20+40 M
90220	AR-7620	NO LONGER AVAILABLE!!
90230	AR-636	VANADIUM PENTOXIDE, 250 G BOTT
90234	AR-608	HALOGEN TRAP Sb METAL
90234	AR-610	HALOGEN TRAP MATERIAL
90235	AR-608	HALOGEN TRAP Sb METAL
90240	AR-263	COPPER ACCELERATOR, 3 POUNDS
90250	AR-508	DISCONTINUED SEE AR2344 AS ALT
90251	AR-639	TIN FLUX 0.5GM PELLETS, 1 POUN
90252	AR-059	Smooth Wall Tin Capsule (H=18m
90256	AR-2183	NICKEL CAPSULE, 6.5 X 10 MM, F
90256	AR-569	NICKEL CAPSULE 4.7 x 9.7, 245m
90256	AR-569-250	NICKEL CAPSULE 4.7 X 9.7, 245M
90258	AR-598	SHORT NIBBLED NICKEL LOW OXYGE
90260	AR-077	IRON CHIP ACCELERATOR, 908 G

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
90270	AR-611	SCHUTZE REAGENT, 100 G PER BOT
90280	AR-076	Tin Metal Accelerator 20 Mesh
90289	AR-170	RARE EARTH COPPER OXIDE, 50 GR
90290	AR-170	RARE EARTH COPPER OXIDE, 50 GR
90291	AR-90291	ACTIVATED CHARCOAL
90300	AR-AEB1300	IRON PHOSPHATE 100G
90301	AR-90300	IRON III PHOSPHATE, 50 G
90330	AR-177	QUARTZ WOOL FINE, 4-12 micron
90331	AR-081	GLASS WOOL, 1 POUND CONTAINER
90332	AR-081	GLASS WOOL, 1 POUND CONTAINER
90339	AR-177	QUARTZ WOOL FINE, 4-12 micron
90340	AR-098	CELLULOSE TRAP MATERIAL 1lb
90341	AR-098	CELLULOSE TRAP MATERIAL 1lb
90341	AR-90341	SEE AR098
90700-1040	AR-091	CALCIUM OXALATE 50g
90710-3010	AR-2092	EDTA, 50 G BOTTLE
90800	AR-073	GRAPHITE POWDER, 100 G
90800	AR-90800	ELTRA STANDARD, GRAPHITE, 100%
90810	AR-1034	CALCIUM CARBONATE STANDARD, 12
90810	AR-1034-100	CALCIUM CARBONATE STANDARD, 12
90810	AR-90810	SEE AR1034
90812-3	AR-90812-3	ELTRA STANDARD, LIMESTONE STAN
90812-3001	AR-4012	C/S IN LIMESTONE STANDARD C= 1
90812-3002	AR-4024	C/S IN LIMESTONE STANDARD C= 1
90815-19	AR-90815-19	CALCIUM CARBONATE 1.95% STD
90816-3001	AR-307	TUNGSTEN CARBIDE C= 6.13% S= 0
90823	AR-90823	DO NOT SELL. CANNOT GET MATERI
90824	AR-90824	SULFANIC ACID STANDARD
90860	AR-211	GLYCINE
90870	AR-5933	COOLING AGENT, 1 PINT CONTAIN
91000-1001	AR-148	COPPER PIN STANDARD, HIGH OXYG
91000-1002	AR-149	COPPER PIN STANDARD, HIGH OXYG
91000-1004	AR-147	COPPER PIN STANDARD, LOW OXYGE
91100-1001	AR-1653	O/N IN STEEL PIN STANDARD O= 0
91100-1002	AR-645	O/N IN STEEL PIN STANDARD O= 0
91100-1003	AR-1651	O/N IN STEEL PIN STANDARD O= 0
91100-1004	AR-1648	O/N IN STEEL PIN STANDARD O= 0
91100-1005	AR-1652	O/N IN STEEL PIN STANDARD O= 0
91100-1006	AR-1656	CURRENT LOT IS LAST OF MATERIA
91100-1007	AR-1650	O/N IN STEEL PIN STANDARD O= 0
91100-1010	AR-1647	O/N IN STEEL PIN STANDARD O= 0
91100-25	AR-91100-25	STEEL PIN STANDARD, 1 GRAM EAC

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
91110-1	AR-91110-1	STEEL STANDARD, HYDROGEN, 1 G
91205-1001	AR-637	H/O/N IN TITANIUM H=0.0013% O=
91205-1003	AR-650	HON IN TITANIUM STANDARD H=0.0
91205-1004	AR-647	H/O/N IN TITANIUM H= 0.0103% O
91205-1005	AR-649	HON IN TITANIUM PIN STANDARD H
91400-1001	AR-546	HYDROGEN IN STEEL PIN STANDARD
91400-1002	AR-555	HYDROGEN IN STEEL PIN STANDARD
91400-1003	AR-556	HYDROGEN IN STEEL PIN STANDARD
91900-1001	AR-4005	ORE STANDARD C= 1.42% S= 1.33%
91900-1002	AR-4006	ORE STANDARD C= 3.75% S= 3.82
92000-1	AR-92000-1	DISCONTINUED
92000-13	AR-92000-13	STEEL CHIP STANDARD FOR CSN, F
92000-14	AR-92000-14	ELTRA CSN STEEL STANDARD, 100G
92000-21	AR-92000-21	ELTRA CSN STEEL STANDARD, 100G
92000-26	AR-92000-26	DISCONTINUED. NO LONGER OFFERE
92000-3	AR-92000-3	ELTRA STEEL CHIP STANDARD
92000-31	AR-92000-31	ELTRA CSH STEEL STANDARD, 100G
92000-35	AR-92000-35	ELTRA CSN STEEL STANDARD, 100G
92000-41	AR-92000-41	ELTRA CSN STEEL STANDARD, 100G
92000-43	AR-92000-43	ELTRA CSN STEEL STANDARD, 100G
92400-3010	AR-673	WROUGHT IRON CHIP STANDARD C=0
92400-3020	AR-946	C/S/N CHIP STANDARD C=0.0050%,
92400-3030	AR-958	C/S/N CHIP STANDARD C= 0.015%
92400-3040	AR-950	C/S/N CHIP STANDARD C= 0.053%
92400-3050	AR-951	C/S/N CHIP STANDARD C= 0.182%
92400-3060	AR-952	C/S/N CHIP STANDARD C= 0.482%
92400-3070	AR-953	C/S/N CHIP STANDARD C= 0.929%
92400-3080	AR-299	CAST IRON C=0.97% S=0.0102%,
92400-3090	AR-301	CAST IRON C=1.99% S=0.0162%,
92400-3100	AR-304	DUCTILE IRON C= 4.21% S= 0.021
92510-30	AR-92510-30	ELTRA COAL STANDARD, 50 GRAMS
92510-45	AR-92510-45	ELTRA COAL STANDARD, 50 GRAMS
92510-67	AR-92510-67	DISCONTINUED
92510-82	AR-92510-82	ELTRA COAL STANDARD, 50 GRAMS
92510-92	AR-92510-92	ELTRA COAL STANDARD, 50 GRAMS
92511-3010	AR-1700	0.35% SULFUR IN COAL STANDARD,
92511-3020	AR-1702	0.73% SULFUR COAL STANDARD
92511-3030	AR-1704	1.09% SULFUR IN COAL STANDARD,
92511-3030	AR-1724	PROX COAL STANDARD, 50 GRAM BO
92511-3040	AR-1705	1.48% SULFUR IN COAL
92511-3040	AR-1706	2.07% SULFUR COAL STANDARD
92511-3040	AR-1726	PROX COAL STANDARD, 50 GRAMS

Cross references: OEM & Alpha Resources Article Numbers

ELTRA	Alpha Resources	Description
92511-3050	AR-1707	2.63% SULFUR IN COAL STANDARD
92511-3060	AR-1709	3.79% SULFUR IN COAL STANDARD
92511-3070	AR-1710	4.25% SULFUR IN COAL STANDARD
92511-3070	AR-1714	1.71% SULFUR IN COAL STANDARD,
92511-3080	AR-1711	5.13% SULFUR IN COAL STANDARD
92511-3080	AR-1712	5.79% SULFUR IN COAL STANDARD
92511-3080	AR-1715	7.06% SULFUR IN COAL STANDARD,
92530	AR-2047	CRUDE OIL, 0.99 WT.%S 100 ML B
92550-3020	AR-2776	ULTIMATE COAL STANDARD, 50 GRA
92550-3040	AR-2781	ULTIMATE ANALYSIS COAL STANDAR
92550-3060	AR-2782	ULTIMATE COAL STANDARD, 50 GRA
92560-3010	AR-2772	ULTIMATE MET COKE STANDARD, 50
92610	AR-241	SILICONE GREASE 150g
95210-67	AR-95210-67	ELTRA COAL STANDARD
95210-82	AR-95210-82	ELTRA COAL STANDARD
AR14016	AR-14016	HOLDER FOR SENSOR AR60352
SVA00040	AR-2783	ULTIMATE COAL STANDARD, 50 GRA
SVA00067	AR-2771	ULTIMATE COKE STANDARD. 50 GRA
SVA00081	AR-1720	PROX COAL STANDARD, 50g
SVA00088	AR-310	CAST IRON C= 3.39% S= 0.080% 1
SVA00092	AR-1622	RESIDUAL OIL ULTIMATE
SVA00093	AR-100	RESIDUAL OIL, C,H,N, BTU 4oz
SVA00097	AR-1683	0.052% S ULTRA LOW SULFUR COAL