

Stearoyl Polyoxylglycerides

Quality Solutions

Application(s): Multiple including ointment base



Authorized Distributor

Stearoyl Polyoxylglycerides are pale yellow, waxy solids. They are dispersible in warm water and in warm paraffin, freely soluble in methylene chloride and soluble in warm methanol. Stearoyl Polyoxylglycerides are used as ointment base, solvents, emulsifying agents and wetting and/or solubilizing agents. USP can support your development and manufacturing activities on stearoyl polyoxylglycerides and medicines containing stearoyl polyoxylglycerides with these existing and upcoming standards.

MONOGRAPH 1

STEAROYL POLYOXYLGLYCERIDES

Status: Official. Official date: MAY 1, 2022 (Recently Revised)

USP Stearoyl Polyoxylglycerides RS

[USP-1621507](#)

08/01/2022. For publications related to **Stearoyl Polyoxylglycerides** in PF online after this date, please visit the integrated [USP-NF/PF online](#).

Disclaimer: USP Quality Solution Sheets are provided as a convenience and for informational purposes only. Labmix24 makes reasonable efforts to provide correct information but assumes no liability for the completeness, timeliness or accuracy of the information contained in each QSS. It is the responsibility of the user to verify the information. The US Pharmacopeia remains the official source of this information.